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**Technical Manual
and
Instructions for Use**

medicap **PRECISE 3000**
homecare GmbH

Electronic gas conserver

Table of Contents

	Page
Contents	Page
1. General explanations and safety instructions	13
1.1. Use	13
1.2. Functional Characteristics	14
1.3. Use- and Safety instructions	14
2. Treatment Preparation	15
3. Running and Settings	15
3.1. Switching on	15
3.2. Power Stage	15
3.3. Select Operation Mode	15
3.3.1. Operation Mode MANUAL	15
3.3.2. Operation Mode AUTOMATIC	15
3.4. Responsiveness	15
3.5. Manual Turn Off	16
3.6. Automatic Turn Off	16
4.. Exchange Bottle	16
5. Maintenance	16
5.1. Cleaning and Disinfection	16
5.2. Check Tightness	16
6. Alarm- and Monitoring Functions	17
7. Technical Data	18
7.1. Electronic Economy Switch	18
7.2. Recommended Accessories	18
8. Guarantee	18

I. General Explanations and Safety Instructions

I.1. Use

The oxygen therapy with the breath-impulse controlled **Precise 3000** is carried out exclusively according to the current medical rules for the use of medical oxygen. By using the **Precise 3000** a multitude of possibilities in the oxygen therapy open up to the user.

- Mobility for the user
- Considerably better use of the oxygen supply, inter alia by automatic adaptation to the O₂ requirement
- Decrease in dessication of the air passages during oxygen inhalation
- Increase of the volume percentage of oxygen in the air
- Best adaptation to the oxygen requirement even during the user's changing performance levels

Latest technical- and manufacturing standards guarantee perfect functioning, high flexibility and utmost user comfort.

Only for patients whose lifefunctions are not directly and permanently dependent on a higher oxygen concentration in the air.

The oxygen inhalation therapy, however, should only be carried out after thorough examination by a physician.

Oxygen for medical purposes is a highly effective drug. Wrong application might cause side effects.

The instructions of the physician are to be strictly observed.

Every physical disorder has to be reported to the physician in attendance immediately.

Contraindications

An oxygen therapy can only be carried out taking special care in the following cases:

- patients of high age
- obesity
- simultaneous ACTH- or glucocorticoid-treatment
- patients with high carbon dioxide concentration in arterial (oxygenrich) blood
- poisoning with substances that reduce breathing.
- breathing control disorders in the central nervous system
- fever

The application of a pure oxygen treatment should not be carried out in case of acute trouble in breathing (respiratory insufficiency on the base of a chronic obstructive emphysema bronchitis) because of the impending decrease in ventilation of the lungs.

Side effects

Taking into consideration the contraindications, side effects under normal oxygen pressure are not to be expected. Patients suffering from insufficient lung ventilation might encounter a rapid increase of the carbon dioxide value inhaling oxygen.

No clinically significant symptoms have been diagnosed during treatments with 50% oxygen for up to seven days. 100% oxygen given for 24 hours, however, leads to cellular and functional lesions of the lungs (cell changes of the respiratory epithelium, secretion inspissation, restriction of ciliar movement, atelectasis as well as changes of the minute volume, carbon dioxide retention and pulmonary vasodilation)

This means that generally during treatments carried out over longer periods of time at 1 atmospheric excess pressure or at even higher oxygen pressures in the air, symptoms of poisoning (hypoventilation, acidosis up to development of an edema of the lungs) are to be expected.

It is to be observed that decreasing the partial pressure too rapidly can lead to a critically insufficient oxygen supply (hypoxemia).

Newborn children receiving highly concentrated (more than 40%)

oxygen treatments for longer periods of time can suffer from lesions to the eyelenses that can cause blindness (retrolenticular fibroplasia).

Apart from that there is danger that bleedings (pulmonial hemorrhage), cellular and/or functional disorders of the lungs (focal atelectasis as well as hyalin membrane lesions with diffuse pulmonary fibrosis) occur.

In order to prevent such a collapse of the lung functions (bronchopulmonial dysplasia), it is imperative to check the oxygen pressure in the arterial (oxygen-rich) blood repeatedly.

1.2. Functional Characteristics

By means of a set pressure reducer the oxygen pressure from the bottle (200 bar) is reduced to a working pressure of 1.6 bar \pm 0.2 bar.

The oxygen is led through a spiral hose from the pressure reducer to the **Precise 3000**. The user is connected to the **Precise 3000** by a nasal canula.

A micro-controller built into the **Precise 3000** regulates the oxygen flow to the user according to two operating modes and sixteen performance steps. The functions of the **Precise 3000** are permanently monitored during use.

The **Precise 3000** only releases the oxygen at the beginning of the inhalation for a certain period of time. Only that oxygen can reach the alveoli and be absorbed by the blood in this way.

Most of the remaining oxygen would be exhaled again.

The **Precise 3000** can be set to the user's individual needs by a medical expert.

The user himself can choose between 16 performance steps according to the respective requirements.

The user must inhale through the nose to make sure the **Precise 3000** can serve properly.

If the device does not receive an inhalation impulse during one minute an acoustic alarm signal is set off and **CHECK NASAL CANULA** appears in the display. If this error is not eliminated the device turns off automatically (see point 3.6).

1.3. Safety Instructions

This technical manual and the instructions for use are part of the unit supplied. It must be available at any time.

Precondition to the agreed use of the **Precise 3000** is thorough knowledge and adherence to these instructions for use.

The instructions given here serve, following the device safety laws, the prevention of dangers by using the device in a not agreed manner and have to be read and observed by all those who use, check and maintain the device.

Do **not** expose the oxygen bottle to any source of heat (heating fan, radiator, oven etc.).

Only qualified service personnel may check the apparatus in case of damages to the plug or mains cable, malfunctioning, general damages from falling down onto the ground or into water.

Humidifiers must not be used.

Observe the important instructions for use.

Only original parts are to be used.

Remove battery when the device is not in use for longer periods of time.

External electromagnetical interference does not endanger the user.

The connections are to be kept absolutely dry and free of grease.

The oxygen bottle is to be protected against falling over. Any damages to the device, the pressure reducer as well as to the oxygen bottle are to be reported to the authorized service personnel.

Keep the device and mountings free of oil and grease (observe the instructions given by the oxygen supplier).

It is absolutely necessary to wash hands before operating the device.

Smoking and open fire generally are strictly prohibited near oxygen leading mountings!

Make sure the oxygen bottle is filled sufficiently.

2. Treatment Preparation

The pressure reducer supplied is to be screwed on manually to the bottle valve with the help of the corrugated swivel nut in clockwise rotations.

Strictly do not use any tools like pliers, wrench or others.

The closed spiral hose is to be connected to the pressure reducer with the help of the swivel nut. By fast connection the spiral hose is to be connected to the connection marked **I** on the Precise 3000.

The nasal canula has to be fixed to the connection marked **↑**.

After this the nasal canula is fitted conveniently. To do this, insert the nasal olives into the nostrils and fix the hose behind the ears with both hands.

3. Operation and Settings

3.1. Switching On

If the apparatus was stored at room temperatures below 10°C the device needs to adapt to the room temperature, otherwise malfunctions might turn up.

Open the valve of the bottle slowly. Switch on the prepared **Precise 3000** by means of the switch **I/O**. There will be a short operational check. The previous settings are indicated and do apply.

The **Precise 3000** can only work correctly if the user inhales through the nose.

3.2. Power Steps

By pressing the buttons marked **< >** the performance steps can be changed in 0.5 l/min steps and perfectly adjusted to the corresponding situation in this way.

The performance steps can be changed during use. See table 1

3.3. Select Operation Mode

Two different operation modes can be set on the Precise 3000. To do this press the two buttons marked **< >** simultaneously and keep them pressed. After about 3 seconds the indication in the display changes from **MANUAL** to **AUTOMATIC** or from **AUTOMATIC** to **MANUAL**.

3.3.1. Operation Mode MANUAL

When set to manual operation mode the amount of oxygen set with the help of the buttons marked **< >** and indicated in the display is led to the user.

The value's setting is fixed and remains, independent from the users' performance levels.

3.3.2. Operation Mode AUTOMATIC

When the breathing frequency is smaller or exactly fifteen breaths per minute the O₂ amount set with the buttons marked **< >** and indicated in the display is led to the user. **X** When the breathing frequency is above 15 and below 30 breaths per minute, the oxygen is increased in a linear way.

● This only happens up the double amount of O₂ set with the buttons marked **< >** not surpassing the maximum amount of 8 l/min!

X The O₂ amount is indicated in the display by full blocks:

currently supplied amount = 2 l/min

● The additional amount supplied is indicated by half blocks:

currently supplied amount = 3 l/min

3.4. Responsiveness

The **Precise 3000** automatically adjusts the sensitivity of the trigger individual to the user.

3.5. Manual Turn Off

The Precise 3000 can be turned off when the treatment has ended pressing the **I/O** button. To do this, press the **I/O** button for approx. 3 seconds. The current settings are stored and an acoustic signal confirms the turn off. Close the oxygen bottle at the bottle valve.

3.6. Automatic Turn Off

When the device is in operation mode and no breathing takes place for one minute **CHECK CANULA** is indicated in the display and an acoustic signal is set off for one minute, interrupted for one minute and set off again for one minute and then the **Precise 3000** turns off automatically storing the current settings. Closing the bottle valve is part of the actions necessary to end the treatment.

4. Changing the Bottle

When the index hand of the pressure manometer falls in the lower range (30bar) of the scale, the bottle needs to be refilled or a new oxygen bottle has to be used in order to ensure perfect functioning of the **Precise 3000**. Never use the oxygen bottle to the very end (see instructions given by the oxygen bottle supplier)

Wash hands thoroughly before each changing of the oxygen bottle. Close the valve before changing.

Switch on the **Precise 3000** pressing **I/O**.

Obtain a pressure free hose system by inhaling a few times. The pressure manometer on the pressure reducer should indicate »0«. Only then the the pressure reducer can be screwd off the bottle valve.

Now change the bottle.

Using a new oxygen bottle the connection to which the pressure reducer is to be screwd onto must be clean. If necessary clean with a dry, grease-free and non-fluffy cloth.

The corrugated swivel nut on the threaded adapter of the pressure reducer is screwd on manual-

ly to the bottle valve turning clockwise.

In any case do not use any tools to do this!

5. Maintenance

5.1. Cleaning and Disinfection

Clean the device now and then with a dry cloth. To clean the mountings (valve, pressure reducer) use a clean, dry cloth only.

The nasal olives and the the nasal canula have to be cleaned after every use inside and outside by wipe-disinfection.

Strictly do not put the entire hose system in disinfecting solution. The user can not guarantee the complete elimination of remaining moisture.

Generally change the nasal canula when changing the user.

5.2. Tightness Check

To check the tightness the device has to be pressure free (see Changing Bottle 4.) After that turn off the Precise 3000. Open the bottle valve slowly. The manometer indicates the respective bottle pressure. Now close the bottle valve. The pressure indicated must not change. Only then the system is tight.

If the pressure indicated is falling there must be a leak in the system. In order to find the leak the screw- and hose connections are daubed off with a soap solution.

If bubbles can be observed either screw fittings again manually or exchange the gasket. Use original parts only.

After that check for tightness again.

Note:

The manufacturer of the pressure reducer demands a general check every 5 years.

Oxygen bottles are subject to TÜV (technical surveillance ass.) checks, i.e. checks in 10-year intervals. The bottles carry, resp. receive a control stamp indicating the date for the next check.

6. Alarm and Monitoring Functions

A built-in microcontroller in the **ECO 300** allows for a permanent monitoring of the most important parameters.

- unhindered oxygen flow
- switching of the magnet valve
- monitoring of breath impulses
- monitoring against hypoventilation
- battery monitoring

When the oxygen flow from the bottle to the **Precise 3000** is impeded an acoustic signal is set off at every inhalation and **NO O₂** is indicated in the display. In that case check the spiral hose connections for correct setting or open the valve, resp. change the bottle in case it is empty.

If the **Precise 3000** is working correctly the switching of the magnet valve at inhalations is indicated by flashing of a * in the upper right hand corner of the display.

When the device is in operation mode and no breathing takes place for one minute **CHECK CANULA** is indicated in the display and an acoustic signal is set off for one minute, interrupted for one minute and set off again for one minute and then the **Precise 3000** turns off automatically storing the current settings. Closing the bottle valve is part of the actions necessary to end the treatment.

When the battery power has become too low **CHANGE BATTERY** appears in the display. In this case change the battery when changing the bottle but not later than five hours after the indication in the display.

If the battery power is too low the **Precise 3000** cannot be used.

In case of **Precise 3000** failure (battery empty, technical defect of the **Precise 3000** etc.) the nasal canula can be connected to the spiral hose with the help of an adapter. A basic alimantation of 2.8 l/min for the user can be ensured in this way.

Operation keys

- I/O** Turn on, turn off the device
< Reduce O₂ flow by 0.5 l/min (minimum 0.5 l/min)
> Increase O₂ flow by 0.5 l/min (maximum 0.5 l/min)

Key combination

- < >** (see 3.3. operation modes)

Indication in the display

Indications in the upper line
 Operation mode **MANUAL** or **AUTOMATIC**
 Switching of the magnetic valve * (in the right corner)

Indications in the lower line

Current amount of O₂

In case of alarm the display changes from the indication of the current amount of O₂ to the alarm message.

Alarm messages:

CHECK CANULA
NO O₂
CHANGE BATTERY

7. Specifications

7.1. Electronic Low Consumption Button

Dimensions: 125 mm x 69 mm x 23 mm without connections

Weight: 220 g without battery

Hose connections: unmistakable, nickel-plated, chromium-plated

Battery: 9 volts-block, type 6LR 61

Working temperature: -5° C to +50° C

Storage temperature: -20° C to +70° C

Classification according MPG: IIa

Triggering: at each breath

Cycle output: 0.5 to 8 l/min correspondingly

Alarms: battery monitoring
missing O₂-supply
no inhalation

7.2 Recommended Equipment

O₂ bottle 2.0 liter

O₂ bottle 0.8 liter

Pressure reducer

Carrier bag for 0.8 l O₂ bottle

Caddy height adjustable + collapsible

Caddy bag for 2,0 l O₂ bottle

Spiral hose

Nasal canula

Adapter for constant flow approx. 2.8 l/min

Battery 6LR61

8. Guarantee

The guarantee is granted for 1 year from the date of delivery. on defects that are due to improper materials or manufacturing faults. Defects covered by the guarantee will be cleared in accordance with our terms and conditions of guarantee service. Medicap does not grant a guarantee if the user endangers the functioning of the device by non-observance of this instruction manual, using the apparatus in a not agreed manner or by third party interference.

Any such cases are the sole liability of the user.

Important

The guarantee is valid only upon presentation of the original invoice/cash receipt.

medicap homecare GmbH

Hoherodskopf Str. 22

D-35327 Ulrichstein

Technical details and illustrations printed in this manual are subject to change.

No copy, translation and reproduction, also in parts, may be made without written permission.

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medicare homecare GmbH