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I. Introduction

I.1. Use

(PRECISE 6000) - for the oxygen enrichment of the air for domestic and clinical oxygen inhalation therapies.

Only for patients whose lifefunctions are not directly and permanently dependent on a higher oxygen concentration in the air.

The oxygen inhalation therapy, however, should only be carried out after thorough examination by a physician.

Oxygen for medical purposes is a highly effective drug. Wrong application might cause side effects.

The instructions of the physician are to be strictly adhered to.

Every physical disorder has to be reported to the physician in attendance immediately.

Contraindications

An oxygen therapy can only be carried out taking special care in the following cases:

- patients of high age
- obesity
- simultaneous ACTH- or glucocorticoid-treatment
- patients with high carbon dioxide concentration in arterial (oxygenrich) blood
- poisoning with substances that reduce breathing.
- breathing control disorders in the central nervous system
- fever

The application of a pure oxygen treatment should not be carried out in case of acute trouble in breathing (respiratory insufficiency on the base of a chronic obstructive emphysema bronchitis) because of the impending decrease in ventilation of the lungs.

Side effects

Taking into consideration the contraindications, side effects under normal oxygen pressure are not to be expected. Patients suffering from insufficient lung ventilation might encounter a rapid increase of the carbon dioxide value inhaling oxygen.

No clinically significant symptoms have been diagnosed during treatments with 50% oxygen for up to seven days. 100% oxygen given for 24 hours, however, leads to cellular and functional lesions of the lungs (cell changes of the respiratory epithelium, secretion inspissation, restriction of ciliar movement, atelectasis as well as changes of the minute volume, carbon dioxide retention and pulmonary vasodilation)

This means that generally during treatments carried out over longer periods of time at 1 atmospheric excess pressure or at even higher oxygen pressures in the air, symptoms of poisoning (hypoventilation, acidosis up to development of an edema of the lungs) are to be expected.

It is to be observed that decreasing the partial pressure too rapidly can lead to a critically insufficient oxygen supply (hypoxemia).

Newborn children receiving highly concentrated (more than 40%)

oxygen treatments for longer periods of time can suffer from lesions to the eyelenses that can cause blindness (retrolenticular fibroplasia).

Apart from that there is danger that bleedings (pulmonary hemorrhage), cellular and/or functional disorders of the lungs (focal atelectasis as well as hyaline membrane lesions with diffuse pulmonary fibrosis) occur.

In order to prevent such a collapse of the lung functions (bronchopulmonary dysplasia), it is imperative to check the oxygen pressure in the arterial (oxygen-rich) blood repeatedly.

I.2. Functional Characteristics

The basic principle of oxygen extraction with the help of oxygen concentrators is the pressure change method following 3 steps:

1st Step: Surrounding air is sucked in by a compressor and is compressed in a container with a molecular sieve

2nd Step: nitrogen molecules are bound at the surface of the molecular sieve and are separated from the oxygen molecules in this way.

The oxygen molecules pass the molecular sieve.

3rd Step: The molecular sieve airs by expelling the nitrogen contained with the help of a part of the extracted oxygen.

By this 3-column method alternately (cyclically) a column is always filled, in the next one the nitrogen is separated from oxygen and in the third one the molecular sieve is aired.

By this method a continuous oxygen flow is provided without an additional oxygen source.

In this way an extremely rapid development of oxygen concentration is attained, right after switching on the device.

1.3. Important Use and Safety Instructions

Precondition for the correct use of the oxygen concentrator (**PRECISE 6000**) is detailed knowledge and adherence of instructions given in this manual.

The apparatus must only be opened by authorized service technicians.

The (**PRECISE 6000**) and all equipment supplied that have contact with oxygen are to be kept free of grease and oil.

Use only original parts.

Do not put the (**PRECISE 6000**) in wet and humid places.

Only qualified service personnel may check the apparatus in case of damages to the plug or mains cable, malfunctioning, general damages from falling down onto the ground or into water.

In case the apparatus cannot be switched off, pull out the plug and call a service technician.

Smoking, open fire and light are not allowed when using the device.

The user is not put at risk by external electromagnetic interferences.

2. Operation

2.1. Preparation

Unfold the carrier handle for convenient transportation.

The (**PRECISE 6000**) should be set up in a well aired, dry, cool and dust-free room.

The (**PRECISE 6000**) should **not** be placed near radiators, fan heaters etc..

The air supply must **not** be obstructed because cooling air is needed for the compressor.

Connect the (**PRECISE 6000**) to the mains cable supplied. Insert plug into mains.

2.2. Switching On

If the apparatus was stored at room temperatures below 10°C the device needs to adapt to the room temperature, otherwise malfunctions might turn up.

Switch the ON/OFF button for automatic check of all functions.

Any malfunctions will be indicated in the display window together with the service telephone number.

2.3. Operate

After switching on the device the flow value administered by your attending physician can be set with the help of the control button of the flow meter (oxygen output in liter/min.).

The middle of the bowl of the flowmeter indicates the set value.

Working hours, operating state and at Option I the treatment hours (ionisation is switched on) are indicated in the display.

When operated with humidifier switch off the hose connection, set the humidifier bottle into the rack and connect the bottle. Then connect the hose to the humidifier.

2.4. Discontinuing Operation

Press the ON/OFF button. The display fades out.

If the device is not used again during the same day disconnect the humidifier bottle. Do not keep the water (risk of germ development).

If not in use over longer periods of time pull the plug from mains.

2.5. Oxygen Concentration Monitoring (OCSI) Option M

A permanent monitoring of the oxygen concentration in the product gas works by using a wear resisting oxygen sensor on an ultrasound basis.

Set values falling below limits are indicated in the display. (see also 4. Troubleshooting Guide).

The sensor does not need new calibration.

2.6. Data Interface Option D

By means of the cable supplied the concentrator can be connected to a computer with SB-interface. This makes the concentrator's service easier.

Apart from that the interface serves as a connection to a low consumption oxygen monitor.

2.7. Ionization Option I

To use the ionization option I unscrew the hose nozzle and connect the ionization supplement in its place. Connect the plug to the relevant bush.

Put the oxygen mask on the ionization head. We recommend to apply an anti-static spray to the mask every time before it is used. Additionally an anti-bacteria filter can be inserted.

Do not use a humidifier during ionization.

The ionization can be activated by pressing I/O on the remote control supplied. The treatment intervals are set with the buttons marked ▽ and △ in 5 minutes steps between 5 and 95 minutes.

The treatment begins by pressing Start. The remaining treatment time is indicated in the display. The treatment can be interrupted by pressing Pause and continued by pressing Start.

An acoustic signal indicates the end of the treatment and End of Treatment appears in the display.

The detachable ionization head should be cleaned once a week or when used by a different patient (see also 5. Cleaning).

3. Diagnosis

An extensive autodiagnosis system has been implemented into the oxygen concentrator (PRECISE 6000).

By pressing the button next to the oxygen outlet on the right the diagnosis is called up.

The displays can be called up by pressing the button.

Display 1: Service Telephone

Display 2: Software Version,
Equipment Version

Display 3: Date, Time

Display 4: Oxygenconcentration at Option M

Display 5: Pressure Evaluation

Display 6: Working Hours Counter

Display 7: Interval Time

Display 8: 1. Error Memory Bank

Display 9: 2. Error Memory Bank

Display 10: 3. Error Memory Bank

4. Troubleshooting Guide

Symptom	Cause	Remedy
Device does not run Display: MAINS Interval Sound	No line voltage Defective fuse	Check if mains cable is connected correctly In case check main fuse Call service
Device is running Display: O2 ← 82%	Flow too high	Set to correct flow
Compressor off Display: O2 ← 72%	Oxygen concentration less 72%	Turn device on and off call service
Device is running Display: O2 Flow too low	Flow regulator turned off Oxygen hose break	Set to correct flow Check hose

Symptom	Cause	Remedy
Compressor off Display: Temperature Interval sound	Device became too hot due to wrong placing Device became too hot because coarse dust filter is blocked Device defect	Let device cool off for 30 min., place correctly, turn on again Clean filter fleece and turn on again Call service
Compressor off Display: Pressure too low Interval sound	Coarse dust filter blocked Intake filter blocked Device defect	Clean filter fleece and turn on again Call service Call service
Compressor off Display: Pressure too high Interval sound	Device defect	Call service
Compressor off Display: Column-valve	Hose connection valve-column not tight	Call service
Device is running Display: Check Head Interval sound	Ionization head not screwd on correctly	Screw on ionization head correctly
Device is running Display: Ions not OK Interval sound	Ionization head is wet Device defect	Dry ionization head Call service

5. Cleaning

Turn off the device and pull the plug from mains before cleaning.

The external box may be cleaned with customary cleaning agents. Do not, however, use any caustic, scratching or acidic cleaning agents.

Do not bring water into the device's interior.

Daily

Rinse the nasal canula with clear, warm water.

When an autoclavable effervescent humidifier is used it should be cleaned, respectively autoclaved, according to the instructions given by the producer of that autoclave effervescent humidify-

er. Non-observance of the instructions can lead to a contamination of the effervescent humidifier.

Option I

Clean Ionization Head

Unscrew ionization head and rinse with clear, warm water. Boil out once a week for 10 to 15 minutes clean in dishwasher at maximum temperature. Leave to dry afterwards. The mask and the ionization head can be treated with the disinfection spray offered. Use only dry ionization heads.

6. Maintenance Intervals

Inspection and maintenance	Once a year or after max. 2800 working hours by the service technician responsible
* Coarse dust filter	Clean filter fleece once a week. If blocked, exchange by the service technician responsible.
* Fine filter	Exchange every six month or after max. 2000 working hours by the service technician responsible.
* Microfilter	Exchange every six month or after max. 2000 working hours by the service technician responsible
Patient bacteria filter	Exchange during inspection and maintenance by the service technician responsible
Inhalation equipment (e.g. nasal canula) and oxygen extension hose	Exchange according to instructions of the attending physician
* Note! All maintenance intervals are guidelines. In different cases these can be shorter due to changing conditions regarding temperature, humidity, dust etc. The filter service life is 20% less in devices with 6 l capacity.	

7. Specifications

	Model (PRECISE 6000) M	Model (PRECISE 6000) M S
Ambient temperature running:	10 to 40° C	10 to 40° C
Dimensions:	55 x 21.5 x 55 cm	55 x 21.5 x 55 cm
Weight:	approx. 21 kg	approx. 21 kg
Line voltage:	230 V +/-10% 50 Hz	230 V +/-10% 50 Hz
Nominal current:	2.0 A	2.0 A
Flow:	0.1 - 5 l/min	0.1 - 6 l/min
Oxygen concentration at 20°C, 1013 hPa 50% rel. humidity	0.1 - 3 l/min 95 +/- 3 Vol.% O2 -4 l/min 90 +/- 3 Vol.% O2 - 5 l/min 82 +/- 3 Vol.% O2	0.1 - 4 l/min 95 +/- 3 Vol.% O2 -5 l/min 88 +/- 3 Vol.% O2 - 6 l/min 80 +/- 3 Vol.% O2
Sound pressure level	38 dB(A)	39 dB(A)
Output pressure	300-350 mbar	300-350 mbar

Option M	
Oxygen measurement	21-98 Vol.% O2 +/- Vol.% O2
Option I	
Ionization capacity	approx. 5 mio. ions/cm ³ /s

8. Equipment supplied

Description	order no.
(PRECISE 6000) M consists of: main unit (oxygencon. monitoring) with mains cable, spare coarse dust filter	006.100
(PRECISE 6000) M I consists of: main unit see PRECISE 6000 M ionization adapter, remote control, tripod	006.110
(PRECISE 6000) M D consists of: main unit see PRECISE 6000 M, interface	006.101
(PRECISE 6000) M S consists of: main unit (higher output, oxygen concentration, monitoring), with mains cable, spare coarse filter	007.100
(PRECISE 6000) M S I consists of: main unit see (PRECISE 6000) M S Ionization adapter, remote control, tripod	007.110
(PRECISE 6000) M S D consists of: main unit see (PRECISE 6000) M S interface	007.101
(PRECISE 6000) M S I D consists of: main unit see PRECISE 6000 M S Ionizationadaptor, remote control, tripod, interface	007.111

Recommended Equipment and Spare Parts:

Nasal canula
Oxygen mask without hose
Oxygen mask with hose
Steril water 500ml Kendall
Fixing device for sterile water 500 ml
Oxygen hose 2 m
Hose connection straight
Angle adaptor with 30 cm silicone hose
Coarse dust filter
Hose nozzle with coupling ring
Power cable
Anti-static spray
Battery filter
Desinfection spray
Tripod
Interface cable Concentrator Computer
Software
Interface Concentrator low consumption monitor

9. Guarantee

The guarantee is granted for 1 year from the date of delivery, on defects that are due to improper materials or manufacturing faults. Defects covered by the guarantee will be cleared in accordance with our terms and conditions of guarantee service. medicap does not grant a guarantee if the user endangers the functioning of the device by non-observance of this instruction manual, using the apparatus in a not agreed manner or by third party interference.

Any such cases are the sole liability of the user.

Important

The guarantee is valid only upon presentation of the original invoice/cash receipt.

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