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(GB) Operation Manual

Précise 6000 Oxygen Concentrator

Picture may differ from original product.

Table of Contents

1. Introduction.....	21
1.1. Use	21
1.2. Functional Characteristics.....	22
1.3. Important Use and Safety Instructions	22
2. Operation.....	24
2.1. Preparation	24
2.2. Switching On	24
2.3. Operate	24
2.4. Discontinuing Operation	24
2.5. Oxygen Concentration Monitoring (OCSI) Option.....	24
2.6. Data Interface Option D.....	25
2.7. Ionization Option I	25
3. Diagnosis.....	26
4. Troubleshooting Guide	27
5. Cleaning	28
6. Maintenance Intervals	28
7. Specifications	29
8. Equipment supplied	30
9. Guarantee.....	31
10. EMC regulations.....	31
11. Oxygen concentrator Précise 6000.....	35
11.1 Interior left side	35
11.2 Interior right side	36
Annex: Declaration of conformity	

1. Introduction

1.1. Use

Précise 6000 - for the oxygen enrichment of the air for domestic and clinical oxygen inhalation therapies. Only for patients whose life functions are not directly and permanently dependent on a higher oxygen concentration in the air. The oxygen inhalation therapy, however, should only be carried out after thorough examination by a physician. Oxygen for medical purposes is a highly effective drug. Wrong application might cause side effects. The instructions of the physician are to be strictly adhered to. Every physical disorder has to be reported to the physician in attendance immediately.

Contraindications

An oxygen therapy can only be carried out taking special care in the following cases:

- patients of high age – obesity
- simultaneous ACTH- or glucocorticoid-treatment
- patients with high carbon dioxide concentration in arterial (oxygenrich) blood
- poisoning with substances that reduce breathing.
- breathing control disorders in the central nervous system
- fever

The application of a pure oxygen treatment should not be carried out in case of acute trouble in breathing (respiratory insufficiency on the base of a chronic obstructive emphysema bronchitis) because of the impending decrease in ventilation of the lungs.

Side effects

Taking into consideration the contraindications, side effects under normal oxygen pressure are not to be expected. Patients suffering from insufficient lung ventilation might encounter a rapid increase of the carbon dioxide value inhaling oxygen.

- No clinically significant symptoms have been diagnosed during treatments with 50% oxygen for up to seven days. 100% oxygen given for 24 hours, however, leads to cellular and functional lesions of the lungs (cell changes of the respiratory epithelium, secretion inspissation, restriction of ciliar movement, atelectasis as well as changes of the minute volume, carbon dioxide retention and pulmonary vasodilation)
- This means that generally during treatments carried out over longer periods of time at 1 atm-superic excess pressure or at even higher oxygen pressures in the air, symptoms of poisoning (hypoventilation, acidosis up to development of an edema of the lungs) are to be expected.
- It is to be observed that decreasing the partial pressure too rapidly can lead to a critically insufficient oxygen supply (hypoxemia).

- Newborn children receiving highly concentrated (more than 40%) oxygen treatments for longer periods of time can suffer from lesions to the eyelenses that can cause blindness (retrolenticular fibroplasia).
- Apart from that there is danger that bleedings (pulmonial hemorrhage), cellular and/or functional disorders of the lungs (focal atelectasis as well as hyalin membrane lesions with diffuse pulmonary fibrosis) occur.
- In order to prevent such a collapse of the lung functions (bronchopulmonial dysplasia), it is imperative to check the oxygen pressure in the arterial (oxygen-rich) blood repeatedly.

1.2. Functional Characteristics

The basic principle of oxygen extraction with the help of oxygen concentrators is the pressure change method following 3 steps:

- **1st Step:** Surrounding air is sucked in by a compressor and is compressed in a container with a molecular sieve
- **2nd Step:** nitrogen molecules are bound at the surface of the molecular sieve and are separated from the oxygen molecules in this way. The oxygen molecules pass the molecular sieve.
- **3rd Step:** The molecular sieve airs by expelling the nitrogen contained with the help of a part of the extracted oxygen.

By this 3-column method alternately (cyclically) a column is always filled, in the next one the nitrogen is separated from oxygen and in the third one the molecular sieve is aired. By this method a continuous oxygen flow is provided without an additional oxygen source.

In this way an extremely rapid development of oxygen concentration is attained, right after switching on the device.

1.3. Important Use and Safety Instructions

Precondition for the correct use of the oxygen concentrator **Précise 6000** is detailed knowlegde and adherance of instructions given in this manual. The apparatus must only be opened by authorized service technicians.

- The **Précise 6000** and all equipment sup-plied that have contact with oxygen are to be kept free of grease and oil. Use only original parts.
- Do not put the **Précise 6000** in wet and humid places.

- Only qualified service personnel may check the apparatus in case of damages to the plug or mains cable, malfunctioning, general damages from falling down onto the ground or into water.
- In case the apparatus cannot be switched off, pull out the plug and call a service technician.
- Smoking, open fire and light are not allowed when using the device.
- The user is not put at risk by external electromagnetic interferences.

Portable and mobile RF communication equipment can affect medical electrical equipment. This may cause a malfunction of the device.

For this reason, radiation sources (devices emitting electromagnetic waves) in the vicinity of medical electrical equipment have to maintain a given minimum distance. See following table.

Radiation sources	Minimum distance from medical electrical equipment, including all connecting cables
Mobile phone (cell phone)	3.3 m
DECT phone (cordless phone)	1.2 m
Bluetooth devices (laptops, cell phone)	0.7 m
Radio remote controls	1.2 m
WLAN devices (e.g.: laptops, repeater, Access Point, Print Server)	2.3 m

The power cord must not exceed a maximum length of 4m.

The use of accessories and cables other than those listed in the instructions for use may lead to increased interference emissions or reduced immunity to interference of the device.

2. Operation

2.1. Preparation

Unfold the carrier handle for convenient transportation. The **Précise 6000** should be set up in a well aired, dry, cool and dust-free room. The **Précise 6000** should **not** be placed near radiators, fan heaters etc.. The air supply must **not** be obstructed because cooling air is needed for the compressor. Connect the **Précise 6000** to the mains cable supplied. Insert plug into mains.

ATTENTION!

If you use Item No. 006.357 FireSafe™ Cannula or Item No. 006.359 FireSafe™ Nozzle please note **Installation Instructions** - FIRESAFE™.

2.2. Switching On

If the apparatus was stored at room temperatures below 10°C the device needs to adapt to the room temperature, otherwise malfunctions might turn up. Switch the ON/OFF button for automatic check of all functions. Any malfunctions will be indicated in the display window together with the service telephone number.

2.3. Operate

After switching on the device the flow value administered by your attending physician can be set with the help of the control button of the flow meter (oxygen output in liter/min.). The middle of the bowl of the flowmeter indicates the set value. Working hours, operating state and at Option I the treatment hours (ionisation is switched on) are indicated in the display. When operated with humidifier switch off the hose connection, set the humidifier bottle into the rack and connect the bottle. Then connect the hose to the humidifier.

2.4. Discontinuing Operation

Press the ON/OFF button. The display fades out. If the device is not used again during the same day disconnect the humidifier bottle. Do not keep the water (risk of germ development). If not in use over longer periods of time pull the plug from mains.

2.5. Oxygen Concentration Monitoring (OCSI) Option

A permanent monitoring of the oxygen concentration in the product gas works by using a wear resisting oxygen sensor on an ultrasound basis. Set values falling below limits are

indicated in the display (see also 4. Troubleshooting Guide). The sensor does not need new calibration.

2.6. Data Interface Option D

By means of the cable supplied the concentrator can be connected to a computer with SB-interface. This makes the concentrator's service easier. Apart from that the interface serves as a connection to a low consumption oxygen monitor.

2.7. Ionization Option I

- To use the ionization option I unscrew the hose nozzle and connect the ionization supplement in its place. Connect the plug to the relevant bush.
- Put the oxygen mask on the ionization head. We recommend to apply an anti-static spray to the mask every time before it is used. Additionally an anti-bacteria filter can be inserted.
- Do not use a humidifier during ionization.
- The ionization can be activated by pressing I/O on the remote control supplied. The treatment intervals are set with the buttons marked and in 5 minutes steps between 5 and 95 minutes.
- The treatment begins by pressing Start. The remaining treatment time is indicated in the display. The treatment can be interrupted by pressing Pause and continued by pressing Start.
- An acoustic signal indicates the end of the treatment and End of Treatment appears in the display.
- The detachable ionization head should be cleaned once a week or when used by a different patient (see also 5. Cleaning).

3. Diagnosis

An extensive autodiagnosis system has been implemented into the oxygen concentrator **Précise 6000**. By pressing the button next to the oxygen outlet on the right the diagnosis is called up. The displays can be called up by pressing the button.

Servicebutton press x times	Designation	Information
1	Language ENGLISH	Select language
2	Softwareversion <<V1.56 –M->>	Software version & equipment e.g. Précise 6000M
3	Serial number	Individual serial number on the equipment
4	O2 Sensor value	Temperature and concentration value
5	Flow Sensor value	Flow value l/min
6	Power switch inverts	Factory adjusted
7	Concentration monitoring: < 80% O2	Factory adjusted
8	Concentration monitoring: < 72% O2	Factory adjusted
9	Battery monitoring	Factory adjusted
10	O2 < 82% Ionization disconnection	Factory adjusted
11	O2 Concentration display	Factory adjusted
12	Flow to small (Monitoring)	Factory adjusted
13	Flow to high (Monitoring)	Factory adjusted
14	Equipment version	Factory adjusted
15	Service telephon number	Servicetelefonnummer des Fachhandelspartners
16	Interval time	Factory adjusted
17	Date	Adjusted
18	Time	Adjusted
19	Working hours counter	1 = Patient 2 = Total
20	min. of pressure values of all 3 columns	Defect of main valve, columns, sensor, compressor etc.
21	max. of pressure values of all 3 columns	Defect of main valve, columns, sensor, compressor etc.
22	Error memory bank of last 10 Alarms	Indicates the error with date and time

4. Troubleshooting Guide

Symptom	Cause	Remedy
Device does not run Display: MAINS Interval Sound	No line voltage Defective fuse	Check if mains cable is connected correctly In case check main fuse Call service
Device is running Display: O2 82%	Flow too high	Set to correct flow
Compressor off Display: O2 72%	Oxygen concentration less 72%	Turn device on and off call service
Device is running Display: O2 Flow too low	Flow regulator turned off Oxygen hose break	Set to correct flow Check hose
Compressor off Display: Temperature Interval sound	Device became too hot due to wrong placing Device became too hot because coarse dust filter is blocked Device defect	Let device cool off for 30 min., place correctly, turn on again Clean filter fleece and turn on again Call service
Compressor off Display: Pressure too low Interval sound	Coarse dust filter blocked Intake filter blocked Device defect	Clean filter fleece and turn on again Call service Call service
Compressor off Display: Pressure too high Interval sound	Device defect	Call service
Compressor off Display: Column-valve	Hose connection valve-column not tight	Call service
Device is running Display: Check Head Interval sound	Ionization head not screwed on correctly	Screw on ionization head correctly
Device is running Display: Ions not OK Interval sound	Ionization head is wet Device defect	Dry ionization head Call service

5. Cleaning

Turn off the device and pull the plug from mains before cleaning. The external box may be cleaned with customary cleaning agents. Do not, however, use any caustic, scratching or acidic cleaning agents. Do not bring water into the device's interior.

Please pay attention to the instructions of the disinfection manufacturer.

For disinfection the following agents can be used:

KORSOLEX FF – Liquid (Bode Chemie, Hamburg)

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Make sure that no Liquid come inside the Précise 6000!

- **Daily**

Rinse the nasal canula with clear, warm water. When an autoclavable effervescent humidifier is used it should be cleaned, respectively autoclaved, according to the instructions given by the producer of that autoclave effervescent humidifier. Non-observance of the instructions can lead to a contamination of the effervescent humidifier.

- **Weekly**

Rinse the nasal canula with clear, warm water (see daily instruction). Afterward to decoct 10 till 15 minutes in water or with dishwasher on high level (**no** fluid or dish liquid).

- **Option I**

Clean Ionization Head Unscrew ionization head and rinse with clear, warm water. Boil out once a week for 10 to 15 minutes clean in dishwasher at maximum temperature. Leave to dry afterwards. The mask and the ionization head can be treated with the disinfection spray offered. Use only dry ionization heads.

6. Maintenance Intervals

Inspection and maintenance	Once a year or after max. 2800 working hours by the service technician responsible
* Coarse dust filter	Clean filter fleece once a week. If blocked, exchange by the service technician responsible
* Fine filter	Exchange every six month or after max. 2000 working hours by the service technician responsible
* Microfilter	Exchange every six month or after max. 2000 working hours by the service technician responsible
Patient bacteria filter	Exchange during inspection and maintenance by the

	service technician responsible
Inhalation equipment (e.g. nasal canula) and oxygen extension hose	Exchange according to instructions of the attending physician

*** Note !**

All maintenance intervals are guidelines. In different cases these can be shorter due to changing conditions regarding temperature, humidity, dust etc. The filter service life is 20% less in devices with 6 l capacity

7. Specifications

	Modell Précise 6000 M	Modell Précise 6000 MS
Ambient temperature running:	10 bis 40 °C	10 bis 40 °C
Dimensions:	55 x 21,5 x 55 cm	55 x 21,5 x 55 cm
Weight:	ca. 21 kg	ca. 21 kg
Line voltage:	230 V+/-10% 50 Hz	230 V+/-10% 50 Hz
Nominal current:	2,0 A	2,0 A
Flow:	0,1 - 5 l/min	0,1 - 6 L/min
Oxygen concentration at 20°C, 1013 hPa 50% rel. humidity :	0,1-3 l/min 95 +/- 3 Vol.% O2 -4 l/min 90 +/- 3 Vol.% O2 -5 l/min 82 +/- 3 Vol.% O2	0,1-4 l/min 95 +/- 3 Vol.%O2 -5 l/min 88 +/- 3 Vol.%O2 -6 l/min 80 +/- 3 Vol.%O2
Sound pressure level:	ca. 38 dB(A)	ca. 39 dB(A)
Output pressure:	300-350 mbar	300-350 mbar
Oxygen measurement:	21-98 Vol.% O2 +/- Vol.-% O2	
Option I		
Ionization capacity:	approx. 5 mio. ions/cm ³ /s	

8. Equipment supplied

Description	order no.
Précise 6000 M consists of: main unit (oxygencon. monitoring) with mains cable, spare coarse dust filter	006.100
Précise 6000 MI consists of: main unit see Précise 6000 M , ionization adapter, remote control, tripod	006.110
Précise 6000 M D consists of: main unit see Précise 6000 M , interface	006.101
Précise 6000 M S consists of: main unit (higher output, oxygen concentration, monitoring), with mains cable, spare coarse filter	007.100
Précise 6000 M S I consists of: main unit see Précise 6000 MS Ionization adapter, remote control, tripod	007.110
Précise 6000 M S D consists of: main unit Précise 6000 MS , interface	007.101
Précise 6000 M S I D consists of: main unit see Précise 6000 MS , Ionization adaptor, remote control, tripod, interface	007.111

Recommended Equipment and Spare Parts:

- Nasal canula
- Oxygen mask without hose
- Oxygen mask with hose
- Steril water 500ml Kendall
- Fixing device for sterile water 500 ml
- Oxygen hose 2 m
- Hose connection straight
- Angle adaptor with 30 cm silicone hose
- Coarse dust filter
- Hose nozzle with coupling ring
- Power cable
- Anti-static spray
- Battery filter
- Desinfection spray
- FireSafe™ Cannula
- FireSafe™ Nozzle
- Tripod
- Interface cable Concentrator Computer
- Software
- Interface Concentrator low consumption monitor

9. Guarantee

The guarantee is granted for defects that are due to improper materials or manufacturing faults. Defects covered by the guarantee will be cleared in accordance with our terms and conditions of guarantee service. medicap homecare does not grant a guarantee if the user endangers the functioning of the device by non-observance of this instruction manual, using the apparatus in a not agreed manner or by third party interference. Any such cases are the sole liability of the user.

Important

The guarantee is valid only upon presentation of the original invoice/cash receipt. Technical details and illustrations printed in this manual are subject to change. No copy, translation and reproduction, also in parts, may be made without written permission.

10. EMC regulations

Regulations and declaration of manufacturer- Electromagnetic compatibility, interference emission		
The Précise 6000 is intended for use in an environment as described below. The operator of the Précise 6000 must make sure that the unit is operated in such an environment.		
Measurement of interference emission	Compatibility	Electromagnetic environment - regulations
HF- emissions acc. to CISPR 11/EN55011	Group1	The Précise 6000 uses HF energy solely for its own function. The HF emission is therefore very low and any interference with adjacent electronic equipment is very unlikely
HF- emissions acc. to CISPR 11/EN 55011	Class B	The Précise 6000 is intended for use in all facilities, including living quarters and other environments that are connected to the public supply network, which also supplies buildings used for the purpose of living.
Emission of harmonics acc. to IEC/EN 61000-3-2	Class A	
Emission of voltage fluctuations / flickers acc. to IEC/EN 61000-3-3	compatible	


Regulations and declaration of manufacturer- Electromagnetic compatibility, interference emission			
The Précise 6000 is intended for use in an environment as described below. The operator of the Précise 6000 must make sure that the unit is operated in such an environment.			
Immunity test	IEC/EN 60601-test level	Compatibility level	Electromagnetic environment - regulations
Discharge of static electricity (ESD) acc. to IEC/EN 61000-4-2	+/- 6kV contact discharge +/- 8kV Air discharge	+/- 6kV contact discharge +/- 8kV Air discharge	Floors should be made of wood or concrete or should be covered with ceramic the floor is covered with synthetic materials, the relative humidity must be at least 30%.
fast transient electric interference / bursts acc. to IEC 61000-4-4	+/- 2kV for mains leads +/- 1kV for input / output leads	+/- 2kV for mains leads +/- 1kV for input / output leads	The quality of the supply voltage should match the typical office or hospital environment.
Surge voltage acc. to IEC/EN 61000-4-5	+/-1kV Push - pull +/-2kV asymmetrical voltage	+/-1kV Push - pull +/-2kV asymmetrical voltage	The quality of the supply voltage should match the typical office or hospital environment.
Voltage dips, shortterm interruptions and fluctuations in supply voltage acc. to IEC 6100-4-11	>5% Ut (>95% dip in Ut) for ½ period 40% Ut (60% dip in Ut) for 5 period 70% Ut (30% dip in Ut) for 25 period >5% Ut (>95% dip in Ut) for 5 s	>5% Ut (>95% dip in Ut) for ½ period 40% Ut (60% dip in Ut) for 5 period 70% Ut (30% dip in Ut) for 25 period >5% Ut (>95% dip in Ut) for 5 s	The quality of the supply voltage should match the typical office or hospital environment.
Magnetic field at a supply frequency (50/60Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at mains frequency must comply with the typical vales found in office or hospital environments.
Note: Ut is the a.c. supply voltage before application of the test levels.			

Regulations and declaration of manufacturer- Electromagnetic compatibility, interference emission			
The Précise 6000 is intended for use in an environment as described below. The operator of the Précise 6000 must make sure that the unit is operated in such an environment.			
Immunity test	IEC/EN 60601- test level	Compatibility level	Electromagnetic environment - regulations
<p>Guided HF-interference according to IEC 61000-4-6</p> <p>Radiated HF-interference in conformance with IEC 61000-4-3</p>	<p>3 V rms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 V rms</p> <p>3 V/m</p>	<p>Portable and mobile radio equipment should not be used closer to the Précise 6000 and the leads, than the safety distance calculated by using the equation applicable for the transmission frequency.</p> <p>Recommended safety distance:</p> <p>$d = 1,2 \cdot \sqrt{P}$</p> <p>$d = 1,2 \cdot \sqrt{P}$; for 80 MHz to 800 MHz</p> <p>$d = 2,3 \cdot \sqrt{P}$; for 800 MHz to 2,5 GHz</p> <p>with P being the rated power of the transmitter in Watt (W) acc. to the specifications of the transmitter manufacturer and d being the recommended safety distance in meters (m)</p> <p>According to an in-situ examination² the field intensity of stationary radio transmitters is with all frquencies lower than the compatibility level³.</p>
Note 1: For 80 MHz and 800 MHz the higher is valid.			
Note 2: These regulations may not apply in all situations. The propagation of electromagnetic waves is influenced by the absorption and reflection by building, objects and persons.			
<p>² The field intensity of stationary transmitters, such as base stations of radiophones and land mobile services, amateur stations, Am- and FM-radio and television stations con theoretically not be exactly predetermined. In order to determine the electromagnetic environment and the effect of stationary HF transmitter, it is recommend to examine the corresponding location. If the detected field intensity at the location of the Précise 6000 exceeds the compatibility level specified above, the Précise 6000 must be examined at a different location, with respect to ist normal operation. If unusual performance characteristics are detected, it may be necessary to apply additional measures, such as reorientation or relocation of the Précise 6000.</p>			
<p>³ Beyond the frequency range from 150 kHz to 80 MHz the field intensity is less than 3 V/m.</p>			

Recommended safety distances between portable and mobile HF communication equipment and the Précise 6000			
The Précise 6000 is intended for operation in an electromagnetic environment with controlled HF interferences. The customer or user of the Précise 6000 can help to avoid electromagnetic interferences by maintaining minimum distances between the portable and mobile HF communication equipment (transmitter) and the Précise 6000 , according to the maximum output power of the communication equipment, as recommended below.			
Safety distance acc. to transmitting frequency m			
Rated power of transmitter W	150 kHz to 80 MHz $d= 1,2 * \text{root from P}$	80 MHz to 800 MHz $d= 1,2 * \text{root from P}$	800 MHz to 2,5 GHz $d= 2,3 * \text{root from P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters for which the rated power is not contained in the table above, the distance can be calculated using the equation for the corresponding column, whereby P represents the rated power of the transmitter in Watt (W), specified by the transmitter manufacturer.			
Note 1			
For calculation of the recommended safety distance for transmitters in the frequency range from 80 MHz to 2,5 GHz an additional of 10/3 was used, in order to reduce the likelihood that a mobile / portable communication unit, that has unintentionally been brought near the patient, will trigger an interference.			
Note2			
These regulations may not apply in all situations, The propagation of electromagnetic waves is influenced by the absorption and reflection by building, objects and persons.			


11. Oxygen concentrator Précise 6000

11.1 Interior left side

Picture may differ from original product.

11.2 Interior right side

Picture may differ from original product.

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TÜV NORD CERT GmbH