**medicap**  
*homecare* ..... GmbH

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**(DE) Gebrauchsanweisung ..... 3**  
**(GB) Instructions for Use ..... 44**

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Edition: Rev.D 2024-10  
Date of issue: 2024-10-08

## Table of contents

1 Introduction .....	46
1.1 Abbreviations.....	46
1.2 How to use the Instructions for Use .....	46
1.3 Significance of warnings and notes.....	46
1.4 Brief description.....	46
1.5 Intended use.....	47
1.6 Functional description.....	47
2. Safety and important information .....	48
2.1 Intended purpose.....	48
2.2 Indications and contraindications .....	48
2.2.1 Side effects and residual risks .....	49
2.3 Intended user and patient target groups.....	49
2.3.1 User responsibility.....	50
2.4 Important instructions for use and safety .....	50
2.4.1 General handling of the device.....	50
2.4.2 Handling recommended accessories .....	51
2.5 Warnings .....	52
2.5.1 Warnings when handling oxygen .....	54
2.6 Expected service life.....	55
2.7 Disclaimer of liability .....	55
2.8 Guarantee and warranty .....	55
2.9 Adresses .....	55
3 Operation .....	56
3.1 Preparation.....	56
3.1.1 Précise 6000 MI / MSI – Connections and control elements.....	56
3.1.2 Installation with Firesafe .....	58
3.2 Switching on the device .....	59
3.3 Diagnostic mode.....	59
3.4 In operation .....	59
3.4.1 Overview of molecular mode therapy .....	60
3.4.2 Overview of ionisation mode therapy.....	62
3.4.3 Setting date & time.....	65
3.4.4 Oxygen therapy - Flow rate .....	65
3.4.5 Night mode .....	66
3.5 Switching off.....	66
3.6 Oxygen Concentration Supervision.....	66
3.7 Service & Diagnosis .....	66
3.7.1 System settings 1: Device information.....	67
3.7.2 System settings 2: Measured values.....	68
3.7.3 System settings 3: Error memory.....	69
4 Malfunctions and their elimination.....	70
5 Cleaning Process.....	73
5.1 General information .....	73
5.2 Cleaning Intervals .....	74
5.3 Cleaning and Disinfection .....	74
6 Maintenance .....	76
7 Pictorial Symbols.....	77
7.1 Definition of the symbols on the device or accessories .....	77
7.2 Software symbols .....	78
8 Technical Data and Specifications .....	80
9 Storage and Scope of delivery.....	81
10 Recommended Accessories, Disposables and Spare Parts.....	82
11 Disposal.....	83

# 1 Introduction

## 1.1 Abbreviations

MD	Medical Device
MDR	Medical Device Regulation (2017/745)
ME	Medical electrical equipment

## 1.2 How to use the Instructions for Use

**Importance of the Instruction for Use** This Instructions for Use is part of the accompanying documents and therefore an essential part of the device. It contains all necessary instructions for operating the device in accordance with Regulation (EU) 2017/745, Annex I: General safety and performance requirements, Chapter III: Requirements regarding the information supplied with the device.

The Instruction for Use must be read thoroughly before the device is put into operation.

The device may only be used by persons who have demonstrably been instructed in its proper use.

**Illustrations** The illustrations used in the documents may differ from the original.

**Changes** Changes to the documents are made as new editions or supplementary sheets. In general, the following applies: Subject to change.

**Reproduction** Reproduction, even in part, is only permitted with written approval.

## 1.3 Significance of warnings and notes

**Warning**

The signal word "Warning" alerts the user to a potentially hazardous situation which, if not avoided, could result in death or serious injury.

**Caution**

The signal word "Caution" draws the user's attention to the fact that in case of failure to observe this warning may result in personal injury, serious injury or damage to the device.

**Note**

The signal word "Note" draws the user's attention to the fact that the following effects are to be expected in case of non-observance:

- Damage to the device may occur.
- Desired functions do not run or do not run correctly.

## 1.4 Brief description

The **Précise 6000** Oxygen concentrator is an electronic device that is used to enrich the breathing air with oxygen during oxygen inhalation therapy at home and in practices. The **Précise 6000 MI / MSI** generates molecular oxygen that can be negatively ionised in the device.

Oxygen for medical purposes is a highly effective medicine. Oxygen inhalation therapy should therefore only be carried out after a thorough medical examination.

The difference between the **Précise 6000 MI** and **Précise MSI** is the maximum flow rate:

**Précise 6000 MI:** 5 l/min

**Précise 6000 MSI:** 6 l/min

## 1.5 Intended use

Precondition for the intended use of the device is the exact knowledge and observance of these Instructions for Use.

The device may only be used in accordance with the general regulations for the installation and operation of medical devices (MPBetreibV – Ordinance on Operators of Medical Devices).

Before using the device, the user must satisfy himself of its functional safety and proper condition.

Read all safety instructions and warnings carefully before using the device. Follow the safety instructions to avoid injuries and life-threatening situations.

When using the accessories, the Instruction for Use of the respective manufacturers and the corresponding warnings must be followed.

If additional devices are connected, the permissible leakage currents may be exceeded. Electromagnetic interference from outside does not pose a risk to the user.

## 1.6 Functional description

Oxygen extraction by the oxygen concentrator works on the principle of molecular sieve technology using the alternating pressure method.

In the

1. Step            ambient air is sucked in by a compressor and compressed in a molecular sieve container (column).
2. Step            nitrogen molecules are bound to the sieve surface and are thus separated from    the oxygen molecules, which can pass through the molecular sieve.
3. Step            the molecular sieve is vented by using some of the recovered oxygen to flush out the nitrogen it contains.

With the 3-column technique used, one column is always filled alternately (cyclically); in another column the nitrogen is separated from the oxygen and in the third column the molecular sieve is flushed. This technology ensures a continuous flow of oxygen without the need of an additional oxygen reservoir. This results in an extremely rapid build-up of oxygen concentration, immediately after switching on the device.

By using a wear-free, ultrasound-based oxygen sensor, the oxygen concentration is permanently monitored.

The negatively charged oxygen ions are generated by means of a corona discharge in the ionisation module. Needle-shaped and current-carrying electrodes release electrons that form negatively charged oxygen ions with oxygen molecules.

The control of the ionisation voltage ensures that no harmful concentrations of ozone (O<sub>3</sub>) and nitrogen oxides can occur.

The ionisation only works when the required high oxygen concentration is present.

## 2. Safety and important information

### 2.1 Intended purpose

The **Précise 6000 MI or MSI** oxygen concentrator is used for oxygen enrichment of the breathing air during oxygen inhalation therapy.

Oxygen inhalation therapy can be used for therapeutic purposes as well as for preventive health care and performance improvement. Therapeutic use of oxygen must take place without restrictions under medical supervision.

Risks were evaluated in the risk management.

**Caution**

The physician's instructions must be followed strictly.

**Caution**

Any physical disorder must be reported to the attending physician immediately.

**Caution**

The effectiveness of the therapy must be checked regularly.

**Warning**

Improper handling and no intended use can cause danger and damage. Therefore, we ask you to read through these Instruction for Use and follow them strictly. Always keep them within easy reach.

### 2.2 Indications and contraindications

#### Indications

Application should only be carried out in accordance with the physician's instructions. The duration and frequency of use should depend on the disease.

Only suitable for patients whose vital functions are not directly and continuously dependent on an increased oxygen concentration in the breathing air.

#### Contraindications

The use of this device is unsuitable in cases of acute respiratory insufficiency (respiratory insufficiency based on chronic emphysema bronchitis) due to the risk of a decrease in lung ventilation.

Contraindications of the accessories can be found in the accompanying documents (e.g. Instructions for Use) of the accessories.

It may only be used on high-risk patients with particular caution and under the special supervision of the attending physician. High-risk patients are persons with:

- Advanced age
- Brain spasm disease
- Immediately after an organ transplant
- Viral infectious diseases
- Pacemaker or cardiac arrhythmia
- Epilepsy
- Somatic diseases
- Obesity
- Simultaneous ACTH or glucocorticoid treatment
- High carbon dioxide concentration in arterial (oxygen-rich) blood
- Poisoning with substances that reduce respiratory activity

- Disorders of respiratory control in the central nervous system
- Fever

### 2.2.1 Side effects and residual risks

If the contraindications are observed, no side effects are to be expected with normal oxygen pressure. Oxygen ventilation of patients with reduced lung ventilation can lead to a rapid increase in carbon dioxide levels.

Side effects may occur as a result of interactions due to improper use of the device, such as the following:

No clinically significant symptoms have been observed during treatments with 50% oxygen for up to seven days. 100% oxygen given over 24 hours, leads to cellular and functional damage of the lungs (cell changes in the alveolar epithelium, thickening of secretion, restriction of ciliary movement, atelectasis and changes in minute volume, carbon dioxide retention and pulmonary vasodilation)

This means that generally during treatments carried out over longer periods of time at 1 at-mosperic excess pressure or at even higher oxygen pressures in the air, symptoms of poisoning (hypoventilation, acidosis up to development of an edema of the lungs) are to be expected.

This means that, as a rule, symptoms of poisoning (hyperventilation, acidosis up to the development of pulmonary oedema) are to be expected after a short period of treatment with 1 atmosphere overpressure over a longer period of time or with even higher oxygen pressures in the breathing air.

It should be noted that reducing the partial pressure too quickly can lead to a life-threatening lack of oxygen (hypoxaemia).

To avoid side effects, medicap homecare GmbH therefore does not recommend treatment for patients at risk.

Newborns who are treated with highly concentrated oxygen (over 40 %) over a longer period of time can suffer damage to the lens of the eye, which can lead to blindness (retrolenticular fibroplasia).

In addition, there is a risk of bleeding (pulmonary haemorrhages), cell and/or functional disorders in the lungs (focal atelectasis) and hyaline membrane damage with diffuse pulmonary fibrosis may occur. To avoid the development of such a collapse in lung function (bronchopulmonary dysplasia), it is essential to repeatedly check the oxygen pressure in the arterial (oxygen-rich) blood during treatment.

For geriatric, paediatric or any other patients who are unable to communicate, discomfort may require additional monitoring. For example, a distributed alarm system to pass on information about the discomfort and/or medical urgency to the responsible nursing staff in order to prevent harm.

### 2.3 Intended user and patient target groups

#### Patient target group

The intended target group are patients with the above-mentioned indications.

Permissible age, gender, duration and frequency of use is to be determined by the attending physician.

The application is not intended for patients with contraindications. These belong to the patient risk groups.

#### **Caution**

The device is not intended for life-supporting or life-sustaining use.

#### **Warning**

Any physical disorder must be reported to the attending physician immediately.

### User group

The product is suitable for home patients and practices.

The Instructions for Use must be studied thoroughly before using the device. The device may only be used by persons who have been instructed in its use.

Patients under the age of 18 or patients with impaired mental or physical condition will only be treated under the supervision of another qualified person.

Patients under 6 years of age may only be treated under the supervision of another qualified person.

Application / measuring in pregnant or breastfeeding women should only be carried out in accordance with the physician's precise instructions

### **2.3.1 User responsibility**

The device may only be operated /used under the operating conditions specified by the manufacturer.

The user may only use the device in the language in which the Instruction for Use is provided.

Note the following when entering parameters:

- The parameters entered must be verified by the user, i.e. the user must check the correctness of the values entered.

The user and/or patient must report all serious incidents related to this device to the manufacturer and/or the competent authority of the Member State in which the user and/or patient is established.

The user does not need any specific qualification, but must follow the requirements for the intended use of the device, depending on the manufacturer's specifications, in these Instruction for Use

The device may only be used by persons who have been instructed in its use.

## **2.4 Important instructions for use and safety**

### **2.4.1 General handling of the device**

**Note**

Precondition for the intended use of the device is the exact knowledge and observance of this Instruction for Use.

---

**Note**

The device must be installed in a well-ventilated, dry, cool and dust-free room.

---

**Note**

The information on the identification plate and the instructions on the labels must be observed.

---

**Hinweis**

If the device has been stored at temperatures below 14 °C, it should be allowed to adjust to room temperature for approx. 60 minutes, otherwise malfunctions may occur.

---

**Note**

Keep the device free from oil and grease.

**Note**

Only original parts should be used. Do not use worn or damaged parts, such as the mains cable.

**Note**

The device must only be opened by an authorised service technician. No changes or modifications may be made to the device.

**Note**

Do not connect the device to the supply voltage via a multiple socket. The mains cable must not exceed a maximum length of 4 m.

**Caution**

The device must be protected from moisture and humidity. DO NOT touch the device when it is wet.

**Caution**

There is no danger to the user from external electromagnetic interference.

**Note**

If possible, the device should not be installed near radiators, fan heaters, etc.

**Note**

Check ALL external components and the packaging for damage. In the case of damage or if the product does not function properly, contact a service technician or the manufacturer so that appropriate repair measures can be carried out.

**Caution**

The air supply must not be impaired, as cooling air is required for the compressor. The minimum distance from cabinets etc. should be 60 cm.

**Caution**

Maintenance and servicing must only be carried out by an authorised service technician.

**Caution**

The oxygen supply setting must be determined individually for each patient with the configuration of the device to be used, including accessories.

**Caution**

If the device is switched on, the fan must run immediately afterwards. If it is not running, please switch the device off again and inform an authorised service technician.

**Caution**

If the device is not used, the **Précise 6000** must be reconnected to the supply voltage every 15 weeks due to the internal battery. After 4 hours, the **Précise 6000** can be disconnected from the supply voltage again. If the device is not used for a longer period of time (more than one week), the mains plug should be disconnected.

**Note**

The device should be used regularly to prevent moisture from forming in the molecular sieve.

## 2.4.2 Handling recommended accessories

**Caution**

The specified accessories, disposables and spare parts have been validated by medicap homecare GmbH for use with the devices.

When using accessories, disposables and spare parts from other manufacturers, accuracy can no longer be guaranteed.

**Caution**

The correct arrangement and positioning of the fixation (e.g.: nasal cannula in the nose or mask on the mouth) is decisive for the amount of oxygen that is delivered to the patient's respiratory system.

**Note**

The accessories, such as the nasal cannula and masks, are "single use" items. This means that they are NOT suitable for repeated use by a single patient or multiple patients.

**Note**

The accessories, such as the nasal cannula and breathing masks, are only to be used for a treatment in accordance with the manufacturer's instructions.

**Note**

The user should be able to hear or feel the flow of gas to the nasal cannula. The hand should be moved back and forth in front of the nasal cannula. If the user does not feel the gas flowing, he must check the connections of the nasal cannula for leaks. Dip the end of the nasal cannula into a cup half-filled with water and check for bubbles.

**Note**

The clinical effect of oxygen > 30 days must be assessed by the attending physician. The possible resulting biological and toxicological risks must be checked and evaluated in advance by the attending physician.

**Note**

When using disposable bubble humidifiers, it must be ensured that they are never operated unopened on the device. Subsequent opening creates negative pressure, which could lead to water ingress and damage to the device. (Flow zero)

## 2.5 Warnings

**Warning**

The device may only be operated by the user under the operating conditions specified by the manufacturer, see chapter 8 Technical data.

**Warning**

Compliance with all Instruction for Use and safety instructions from Chapter 2 is mandatory.

**Warning**

Smoking, open light and fire are prohibited during use.

**Warning**

The device must not be used if:

- the mains plug or cable is damaged
- functions run incorrectly
- the device has been dropped
- the device is damaged
- liquid has penetrated the device

In these cases, the device must be checked by an authorized service technician and repaired if necessary.

**Warning**

The device must only be opened by an authorised service technician.

**Warning**

Multiple sockets or extension cables must not be used.

**Warning**

The device should not be placed near radiators, fan heaters, etc.

**Warning**

The use of accessories and cables other than those listed in the Instruction for Use can lead to increased interference emission, reduced interference immunity of the device or impairment of the application.

**Warning**

DO NOT use the device in an environment where contamination, smoke, gases, flammable anaesthetics, cleaning agents or chemical vapours are present.

**Warning**

Close supervision is required when this device is used in the vicinity of children or disabled persons.

**Warning**

Portable and mobile HF communication equipment can affect medical electrical devices. This can lead to malfunctions in the devices.

For this reason, radiation sources (devices that emit electromagnetic waves) must be kept at a certain minimum distance from operating medical electrical devices. See table below:

**Radiation sources**

**Minimum distance from a medical device including all connection cables**

Mobile phone	3.3 m
DECT telephone (cordless telephone)	1.2 m
Bluetooth devices (laptops, mobile phones)	0.7 m
Remote controls	1.2 m
WLAN devices (e.g.: laptops, repeater, access point, print server)	2.3 m

**Warning**

Patients who cannot see or hear alarm signals or cannot communicate their discomfort must be supervised.

**Warning**

A change in altitude above sea level can affect the performance of the device. Consult your physician or contact the manufacturer.

**Warning**

Stow and lay electrical cables properly to avoid potential tripping hazards.

**Warning**

Risk of infection if the device is used again!  
If the device is used by several patients, infections can be transmitted to the next patient.

When the device is used again, it must be cleaned and disinfected, see chapter 5.

**Warning**

The ME equipment has no defibrillation protection.

**Warning**

Modification of the ME equipment is not permitted.

**Warning**

The device may not be modified without the manufacturer's permission.

**Warning**

If the ME equipment is modified, suitable inspections and tests must be carried out to ensure continued safe use, suitable inspections and tests must be carried out.

## 2.5.1 Warnings when handling oxygen

### Warning

There is a fire risk associated with the enrichment of oxygen during oxygen therapy. Do not use the oxygen concentrator or accessories near sparks or open flames.

### Warning

To ensure that the therapeutic amount of oxygen delivery is delivered according to your medical condition, the **Précise 6000 MI / MSI**:

- should only be used with settings that have been individually determined or prescribed for you and your current activity level with your accessories;
- may only be used with the specific combination of parts and accessories that meet the specification of the concentrator or accessory manufacturer.

### Warning

Before and during oxygen therapy, use only water-based lotions or ointments that are compatible with oxygen. Do not use petroleum or oil-based lotions or ointments to avoid the risk of fire or burns.

### Warning

Fittings, connections, tubings or other accessories of the oxygen concentrator must not be lubricated in order to avoid the risk of fire or burns.

### Warning

To ensure proper function and to avoid the risk of fire or burn injuries, spare parts recommended by the manufacturer must be used.

### Warning

Oxygen facilitates the sparking and spread of fire. Do not leave the nasal cannula or mask on bed covers or chair cushions when the oxygen concentrator is switched on and not in use; The oxygen increases the flammability of the material and makes it flammable. Switch off the oxygen concentrator when it is not in use to prevent oxygen enrichment.

### Warning

If this device is used above an altitude of 2000 meters above sea level, or outside a temperature of 15°C - 40°C or at a relative humidity of more than 65%, a negative influence on the volume flow and the percentage of oxygen and consequently an impairment of the quality of the therapy is to be expected.

### Warning

If you experience discomfort or a medical emergency during oxygen therapy, seek medical attention immediately to avoid harm.

### Warning

For geriatric, pediatric or any other patients who are unable to communicate, complaints may require additional monitoring and/or a distributed alert system to relay information about the complaint and/or medical urgency to the appropriate caregiver to prevent harm.

### Warning

Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking or open flames in the room where the oxygen concentrator or oxygen-containing accessories are located. If you want to smoke, you must always switch off the oxygen concentrator, remove the nasal cannula and leave the room in which either the nasal cannula or the mask or the oxygen concentrator is located. Unless you are able to leave the room, you must wait 10 minutes after switching off the oxygen concentrator before smoking.

### Warning

Open flames during oxygen therapy are dangerous and may lead to fire or death. Open flames within 2 m of the oxygen concentrator or oxygen-carrying accessories are not permitted.

### Warning

Oxygen tubings, cables and the concentrator must not be covered by blankets, covers, seat cushions, clothing, etc. or come close to heated or hot surfaces, such as radiators and similar electrical devices.

## 2.6 Expected service life

If maintenance is carried out to the prescribed extent and at the specified intervals, the device will operate safely during the intervening period.

With each maintenance and Technical Safety Control, the “expected service life” according to IEC 60601-1 will therefore be prolonged until the next prescribed Technical Safety Check.

## 2.7 Disclaimer of liability

The device has been approved for use with certain accessories and spare parts (see Chapter 10).

The manufacturer accepts no responsibility or liability for personal injury or other damage and excludes any warranty for damage to the device resulting from the use of non-approved or unsuitable accessories or spare parts

## 2.8 Guarantee and warranty

medicap homecare GmbH provides a two-year guarantee from the date of delivery for defects caused by material or manufacturing faults.

Defects covered by the warranty will be rectified in accordance with the warranty conditions.

In addition, medicap homecare GmbH does not provide a guarantee if the operator jeopardizes the functions of the device by failing to observe these Instructions for Use, by using the device for purposes other than those for which it is intended or by tampering with it.

medicap homecare GmbH does not automatically consider the owner of the device to be authorised to carry out service work.

### Note

The guarantee can only be claimed in conjunction with the proof of purchase.

## 2.9 Adresses

Manufacturer:

medicap homecare GmbH  
 Hoherodskopfstr. 22  
 35327 Ulrichstein  
 Germany  
 Tel.: +49 6645 970-0  
 Fax: 06645/970-200

## 3 Operation

### 3.1 Preparation

Please check that all the required elements are present before assembly. Keep the provided packaging safe and dry for the envisaged servicing.

Unfold the integrated carrying handle for manoeuvring the **Précise 6000 MI / MSI**. The **Précise 6000 MI / MSI** must be installed in a well-ventilated, dry and dust-free room. The **Précise 6000 MI / MSI** must not be installed near radiators, fan heaters etc. The air supply should not be impaired as cooling air is required by the compressor.

The **Précise 6000 MI / MSI** must be connected using the power cable, as provided. Insert the mains plug into a socket outlet, do not connect to the supply voltage via a multiple socket outlet.


#### Caution

If the mains voltage is interrupted or the fuse is defective, the ON/OFF button lights up three times and a signal sounds.


#### 3.1.1 Précise 6000 MI / MSI – Connections and control elements

Front view



Rear view:

Detail view mains socket:

Detail view connections:

### 3.1.2 Installation with Firesafe

Screw the Firesafe™ valve adapter hand-tight to the **Précise 6000 MI / MSI** Oxygen outlet.

Push the ionisation module from below into the holder (black knob) on the oxygen concentrator and with the aluminium connector onto the Firesafe™ valve.

Insert the ionisation module plug into the device.

### 3.2 Switching on the device

If the device has been stored or transported at temperatures below 14 °C, the device must first acclimatise to room temperature, otherwise malfunctions may occur.

To switch on, press the ON/OFF button; a dual acoustic signal sounds, the device starts and carries out a self-test (diagnostic mode); this takes approx. 45 seconds. In the event of a fault, an audible alarm sounds and the error code is shown on the display.


In the event of a serious fault, the device automatically switches to error mode or a shutdown is triggered.

#### Caution

In diagnostic mode, the flow value icons are inactive and cannot be changed.

### 3.3 Diagnostic mode

### 3.4 In operation

The molecular oxygen supply starts after the diagnostic mode automatically.

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Molecular mode therapy is automatically activated on delivery.

#### Note

The timer  The timer is only used to time the application.

Generally, the use of disposable bubble humidifiers with sterile water is recommended for molecular oxygen inhalation. For reusable bubble humidifiers, the humidifier bottles must be cleaned in accordance with the manufacturer's instructions.

The desired molecular mode Standard or ionisation mode Standard can only be preset by a service technician:

**Note**

The desired molecular mode Standard or ionisation mode Standard can only be preset by a service technician:

- Standard: Oxygen supply without timer function. An application duration cannot be set.

**Note**


Molekular mode therapy or ionisation mode therapy:  
Application duration can be set with an acoustic signal after the set time has elapsed.

### 3.4.1 Overview of molecular mode therapy

The application duration can only be set before starting the application; once the application has been successfully started, the application duration can no longer be set.

- Maximum 2 hours
- Minimum 5 minutes
- Adjustments possible at 5-minute intervals
- 

The user can pause and stop the application using the symbols  and .

After the set duration of use, an acoustic signal sounds and the therapy ends. The molecular oxygen supply remains active.


### 3.4.2 Overview of ionisation mode therapy

The application duration can only be set before starting the application; once the application has been successfully started, the application duration can no longer be set.

- Maximum 2 hours
- Minimum 5 minutes
- Adjustments possible at 5-minute intervals

The user can pause and stop the application using the symbols  and .

Molecular mode is always active until ionisation is activated. Ionisation is activated by pressing the “ Start application” symbol.

**Note**

If the ionisation module is correctly connected, the active ionisation mode  automatically appears on the display. The active molecular mode symbol  appears if no ionisation module is connected or if it is not connected correctly.

After the set application time, an acoustic signal sounds and the ionisation therapy ends.

The therapy can also be ended by pressing the “ Stop application” symbol.

Once ionisation therapy has ended, the molecular oxygen supply remains active. This means that the molecular oxygen supply continues.

### Overview of status notifications


Application running	 <b>Molecular running</b>
Application paused	 <b>Molecular paused</b>
Molecular“ application successfully completed	 <b>Molecular completed</b>
„Ionisation“ application successfully completed	 <b>Ionisation completed</b>

### 3.4.3 Setting date & time

To access the "Date & Time" menu, tap on the area highlighted in red. The following menu guide appears:

Use the  and  to set the date and time. The setting once effected needs to be confirmed by pressing  and stored. After successful setting, this must be confirmed with 0 and saved. The system software is then re-booted.

By pressing  you can quit the menu without any changes; no re-booting of the system software is needed.

#### Caution

The date and time can only be set if no application is running!


### 3.4.4 Oxygen therapy - Flow rate

After switching on and successfully passing the diagnostic mode the SET flow rate recommended by your physician (oxygen output in litres/min) can be set by pressing the  and  symbols.

The change is made in 0.2 l/min. steps; the required SET flow rate is shown for approx. 5 seconds and then automatically changes to the ACTUAL flow rate.

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### 3.4.5 Night mode

Night mode can be set on the device. The time after which the device switches to night mode can be set between 1 and 10 minutes. The time setting must be confirmed with .

If night mode is active, the display no longer lights up. Night mode can be deactivated by simply touching the display.

### 3.5 Switching off


Press ON/OFF button for approx. 3 seconds. The display switches off and an acoustic signal sounds. If there is no more operation on the same day, the humidifier bottle used must be cleaned. Do not store the water (germ formation). If the device is not to be used for a longer period of time, the plug must be removed from the socket.

### 3.6 Oxygen Concentration Supervision

The oxygen concentration in the product gas is permanently monitored using a wear-free, ultrasound-based oxygen sensor. If limit values are not reached, these are shown on the display. The sensor does not need to be recalibrated.


### 3.7 Service & Diagnosis

The **Précise 6000 MI / MSI** oxygen concentrator features comprehensive self-diagnostics.

To access the "System settings", tap on the area marked in red. After exiting the system settings automatically or manually, the diagnostic mode is executed. After approx. 30 seconds of inactivity, the software automatically exits the system settings.

### 3.7.1 System settings 1: Device information


**Device information**

Version software	1.10 (12.09.2022)
Version Display-Controller	Apr 07 2022
Version power controller	Mrz 04 2022
Version DB-Text	v1.6 (26.03.2024)
Version DB-Data	v1.2 (13.07.2021)
Software language	
Service telephone	
Device ID	600.100.00.0000
Compressor no.	062002000847
Operating hours meter	00217:24
IP-address	
Night mode (min)	

Annotations:

- Volume On / Off (points to the speaker icon)
- Language selection (points to the language dropdown menu)
- Overview 2 - Measurements (points to the information icon)

If the language has been changed, the following message is displayed:


### 3.7.2 System settings 2: Measured values

Measured values 			
	Column 1	Column 2	Column 3
<b>Min</b>	0.72 bar	0.72 bar	0.73 bar
<b>Max</b>	1.04 bar	1.06 bar	1.06 bar
<b>Interval time</b>		3.7 sek	
<b>O2 concentration</b>		91.9 Vol%	
<b>Temperature</b>		22.7 °C	
<b>Flow</b>		4.0 l/min	
<b>Rechargeable battery voltage</b>		14,13 V	

Overview 3 -  
Error memory

### 3.7.3 System settings 3: Error memory

Exiting the system setting

## 4 Malfunctions and their elimination

Failure symptom	Cause of error	Remedy
 <p>Service telephone</p> <p><b>Mains voltage interrupted!</b></p>	<p>No power or device fuse defective</p> <p>Device switches off automatically</p>	<p>Check whether the power cable is connected correctly, check the device fuse or main fuse, see picture „Detail view mains socket“</p> <p>Call Service</p>
 <p><b>Oxygen concentration &lt; 82%</b></p>	<p>Oxygen concentration &lt; 82%</p>	<p>Reduce flow</p> <p>Call Service if the message is repeated</p>
 <p>Service telephone</p> <p><b>Oxygen concentration &lt; 72%</b></p>	<p>Oxygen concentration &lt; 72%</p> <p>Device automatically switches to error mode</p>	<p>Switch device off and on again</p> <p>Call Service if the message is repeated</p>
 <p><b>Oxygen flow too low</b></p>	<p>Oxygen tube to the mask kinked</p>	<p>Check the tube or mask</p>

 <p>Service telephone</p> <p><b>Temperature too high!</b></p>	<p>Device has overheated 1) Device was positioned in the wrong place, 2) Coarse dust filter is clogged</p> <p>Device switches automatically to error mode</p>	<p>Allow device to cool down for 30 min, set up correctly and switch on again, clean the filter fleece and switch on the device again</p> <p>Call Service if the message is repeated</p>
 <p>Service telephone</p> <p><b>Pressure too low!</b></p>	<p>1) Coarse dust filter is clogged 2) Suction filter is clogged 3) Device defective</p> <p>Device automatically switches to error mode</p>	<p>Clean the filter fleece and switch it on again</p> <p>Call Service if the message is repeated</p>
 <p>Service telephone</p> <p><b>Pressure too high!</b></p>	<p>Device defective (column pack too damp)</p> <p>Device automatically switches to error mode</p>	<p>Call Service</p>
 <p>Service telephone</p> <p><b>Column / valve</b></p>	<p>Tube connection valve - column leaking</p> <p>Device automatically switches to error mode</p>	<p>Call Service</p>
 <p><b>Check ionisation head!</b></p>	<p>Ionising head not screwed on correctly</p>	<p>Check ionisation head, screw on correctly and hand-tight</p>

 <p>Ionisation voltage too low</p>	<p>Ionisation voltage too low</p>	<p>Ionising head is too damp Dry ionising head</p>
 <p>Ionisation finished! Oxygen concentration &lt; 82%</p>	<p>Oxygen concentration &lt; 82%</p>	<p>Reduce flow Call Service if the message is repeated</p>

## 5 Cleaning Process

Equipment cleaning is a process that a used product undergoes so that it can be safely reused. This procedure includes cleaning and disinfection.

In addition, a test and restoration of the technical and functional safety of the product used must be carried out.

### Warning

Risk of infection if the device is used again! If the device is used by several patients, infections can be transmitted to the next patient. If the device is used again, it must be cleaned and disinfected.

### 5.1 General information

### Note

Wear suitable protective equipment during cleaning / disinfection.

### Caution

Switch off the device and disconnect it from the power supply before starting cleaning.

### Caution

The cleaning and disinfecting agent must be used in accordance with the manufacturer's instructions, see the manufacturer's Instructions for Use.

### Caution

Do not allow any liquid to enter the interior of the device.

### Caution

DO NOT allow detergent to drip into the air inlet or outlet ports.

### Caution

After cleaning, disinfection and a functional check, the device is suitable for reuse with further patients.

### Warning

Risk of injury due to electric shock! Penetrating liquids from cleaning agents and disinfectants can cause a short circuit, injure the user and damage the device.

### Caution

Ensure spatial separation of equipment that is still contaminated and equipment that has already been disinfected (such as work surfaces and storage rooms, etc.).

### Caution

Do not open the device! Repairs may only be carried out by qualified technicians.

### Caution

Do not pour liquid over the device and components or immerse them in liquids.

### Caution

DO NOT clean the housing with concentrated chlorine-based products (ethylene chloride), oil-based products (Pine-Sol®, Lestoil®), or other aggressive chemical cleaning agents. Use commercially available cleaning agents to clean the outer housing, but do not use abrasive, scratching or corrosive agents.

### Caution

No lubricants other than those recommended by medicap homecare GmbH may be used.

## 5.2 Cleaning Intervals

Intervals	Action
Daily	- Reusable (bubble) humidifier, in accordance with the manufacturer's instructions
Weekly	- Coarse dust filter - Ionisation head
Monthly	- Clean device
Annual	- Replace fine filter - Replace the microfilter
If required	- Clean device - Ionisation head - Replace / exchange filter / oxygen connection tube - Replace / exchange accessories
For patient change	- Before using the device again, have it cleaned and disinfected by a specialised dealer. - After every course of treatment and every time the patient is changed, replace the oxygen connection tube. - Before using the ionising head again, have it cleaned and disinfected by a specialist dealer. - Replacing the bacteria filter disc.

## 5.3 Cleaning and Disinfection

Component / part	Cleaning (manual)	Disinfection (manual)
Housing incl. handle	Wipe with a damp cloth; Use water or commercially available cleaning agents, but do not use abrasive, scratching or corrosive agents, Contact time max. 2 min, allow to dry	Wipe disinfection (see recommended disinfectants below), allow to dry.
ON /OFF Button or other control elements	Wipe with a damp cloth; Use water or commercially available cleaning agents, Contact time max. 2 min, allow to dry	
O <sub>2</sub> Outlet	Clean the outlet with an alcohol swab and air dry	Wipe disinfection (see recommended disinfectants below), allow to dry.
Display (Touchscreen)	Wipe with a damp cloth; Use water or commercially available cleaning agents, Contact time max. 2 min, allow to dry	Wipe disinfection (see recommended disinfectants below), allow to dry.
Power supply unit / power cable		

Ionisation module	Wipe with a damp cloth; Use water or commercially available cleaning agents, Contact time max. 2 min, allow to dry	Wipe disinfection (see recommended disinfectants below), allow to dry.
Bracket		
Ionisation head	Unscrew the upper part of the ionisation head and remove the internal bacteria filter disc. The upper part must be cleaned with clear warm water and can be disinfected according to the manufacturer's instructions. The ionisation head must be dried. A new bacteria filter disc must be inserted when the unit is used again.	Disinfection according to the manufacturer's instructions (see disinfectant recommendation below).  Allow to dry.
Reusable bubble humidifiers	For reusable bubble humidifiers, the humidifier bottle must be cleaned in accordance with the manufacturer's instructions	For reusable bubble humidifiers, the humidifier bottle must be disinfected in accordance with the manufacturer's instructions
Nasal cannula / mask	Single-use items are to be disposed of, no cleaning	Single-use items are to be disposed of, no disinfection
Coarse dust filter	Rinse with water and be sure to dry! Replace if clogged	No disinfection
Bacteria filter disc	No cleaning - Replacement	No disinfection

Cleaning of the devices is always done manually. The cleaning agent is to be applied on a cloth (cloths with low particle release are to be used).

Brushes, spray guns or similar are not recommended.

#### Recommendation for cleaning agents and disinfectants

The following disinfectants have been validated by medicap homecare GmbH and can be used:

- Biguanid Fläche N (Dr. Schumacher)
- Descogen Liquid
- microzid AF wipes (Schülke)
- 5% Helipur H plus N from-B Braun for the ionising head

Any other disinfectant used must be listed in the VAH list and be compatible with the materials to be treated.

Drying is carried out in fresh air without certain drying media (such as hair dryers). Wiping with cloths (for drying) is not recommended.

## 6 Maintenance

**Caution**

All maintenance intervals are guide values. In individual cases, these may be shorter due to changes in environmental conditions with regard to temperature, humidity, dust, etc.

---

**Caution**

For devices with a capacity of 6 l/min, the filter service life is reduced by up to 20 %.

---

**Caution**

Maintenance may only be carried out by an authorised service technician.

---

**Caution**

During maintenance, only monitored test equipment and measuring devices suitable for the intended purpose may be used.

The device and the sensor must be checked and serviced annually by the supplier or manufacturer or by authorised trained personnel.

For commercially used devices, a safety check must be carried out.

The following must be replaced during annual maintenance by the responsible service technician:

Fine filter, microfilter, patient bacteria filter.

The device must not be maintained or serviced during use on the patient!

---


**Caution**

Works on the opened device may only be carried out by an authorised service technician.

---

## 7 Pictorial Symbols

### 7.1 Definition of the symbols on the device or accessories

	Attention!
	No smoking
	No open flames
	Follow the Instructions for Use
	Follow the Instructions for Use
	Fragile
	Keep away from sunlight
	Store in a dry place
	Protection class II
	Single-use; Do not reuse
	Contain Latex
	Batch number
	Serial number
	The device must not be disposed of with household waste
	CE- Mark 0044; Notified body: TÜV NORD CERT GmbH 0044
	Degree of protection against electric shock: Type BF

	Alternating current
	Dangerous electric-voltages
	Manufacturer
	Date of manufacture
	Medical device
	UDI number - Unique device identifier
	Article number

## 7.2 Software symbols

	Duration time		Exit / Back
	Increase		Reduction
	Start application		Pause application
	Stop application		Ionisation mode
	Setting figure – upwards		Setting figure – downwards
	Active		Inactive
	Accept		Cancel
	Molecular mode		System settings
	Warnings		Flow rate

 <p>Sound An / Sound on</p>	 <p>Sound off</p>
 <p>Increasing the illuminance</p>	 <p>Illuminance reduction</p>
 <p>Fault notification - Device switches to the error mode</p>	 <p>Mode – ACTIVE</p>

## 8 Technical Data and Specifications

Functionality is guaranteed at room temperature, normal humidity (10 - 65 %), good lighting conditions (at least 750 lux) and no background noise.

The ambient volume should not exceed room volume, up to a maximum of 55 dB.

All the parameters required for the device to function properly are entered in the following table:

	<b>Modell Précise 6000 MI</b>	<b>Modell Précise 6000 MSI</b>
MDD 93/42 EEC classification:	Ila	Ila
Degree of protection against penetration of liquids:	IP21 Drip water-proof	IP21 Drip water-proof
Operating temperature:	15 bis 40 °C	15 bis 40 °C
Dimensions:	55 x 21,5 x 55 cm	55 x 21,5 x 55 cm
Weight:	ca. 21 kg	ca. 21 kg
Operating voltage:	230 V+/-10% 50 Hz	230 V+/-10% 50 Hz
Rated current:	2,0 A	2,0 A
Power consumption:	350 W	350 W
Flow:	1 - 5 l/min	1 - 6 l/min
Oxygen concentration:	1-3 l/min 95 +/- 3 Vol.% O <sub>2</sub> -4 l/min 90 +/- 3 Vol.% O <sub>2</sub> -5 l/min 82 +/- 3 Vol.% O <sub>2</sub>	1-4 l/min 95 +/- 3 Vol.% O <sub>2</sub> -5 l/min 88 +/- 3 Vol.% O <sub>2</sub> -6 l/min 80 +/- 3 Vol.% O <sub>2</sub>
Sound level (device):	ca. 39 dB(A)	ca. 39 dB(A)
Outlet pressure:	max. 350 mbar	max. 350 mbar
Oxygen measurement:	21-96 Vol.% O <sub>2</sub> +/- 3 Vol.-% O <sub>2</sub>	21-96 Vol.% O <sub>2</sub> +/- 3 Vol.-% O <sub>2</sub>
Operating heighth:	0 bis 2000 m	0 bis 2000 m
Fuse:	2 x T2,5 A	2 x T2,5 A

## 9 Storage and Scope of delivery

The device is delivered safely packed in the box. It is transported and stored in the original box. Store the device, before and after use, under the prescribed environmental conditions.

If the device has been used, the device should be cleaned and disinfected before storage.

Picture	Description	Article number	Quantity
	<b>Précise 6000 MI</b> Oxygen concentrator	600.110	1
	<b>Précise 6000 MSI</b> Oxygen concentrator	700.110	1

Article number	Description
600.110	<b>Précise 6000 MI</b> consisting of: Oxygen concentrator with Instructions for Use, Power cable 4m, Ionisation module external, Bacteria filter discs (5 pcs.), FireSafe™ Nozzle with aluminium adapter, Replacement coarse dust filter, 2x wheel for housing, Allen key for wheels 4mm, Antistatic spray, Mask for ionisation, Grab rail ION
700.110	<b>Précise 6000 MSI</b> consisting of: Oxygen concentrator with Instructions for Use, Power cable 4m, Ionisation module external, Bacteria filter discs (5 pcs.), FireSafe™ Nozzle with aluminium adapter, Replacement coarse dust filter, 2x wheel for housing, Allen key for wheels 4mm, Antistatic spray, Mask for ionisation, Grab rail ION

## 10 Recommended Accessories, Disposables and Spare Parts

The accessories, disposables or spare parts listed below have been validated by medicap homecare GmbH for use with the Précise 6000 medical device and are recommended.

**Caution**

The accessories such as the nasal cannula or masks are "single use" items. They are NOT suitable for repeated use by a single patient as well as by multiple patients in the physician's office/clinic.

**Caution**

The use of accessories, spare parts or disposables that have NOT been approved may result in impaired performance.

**Caution**

The use of the "FireSafe" accessory protects the patient during the oxygen supply in the event of a fire by preventing the spread of flames.

Article number	Article Description
006.900	Nasal cannula for adult, straight prongs
006.901	Oxygen respiratory mask + 2,1 m extension tube
006.907	Bubble humidifier, disposable 340ml
006.359	FireSafe Nozzle with aluminium adapter
006.322	O <sub>2</sub> -outlet tube connector with cap nut
6079	Straight connector Norma GS 5 plastic
006.905	Extension tube "O <sub>2</sub> " 2.1 m
006.906	Extension tube 10 m
006.183	Coarse dust filter
006.184	Fine filter
006.185	Micro filter, Filter ISO Gard; HEPA light
006.902	Oxygen respiratory mask for Ionisation
016.025	Antistatic spray for Ioniser, 250 ml ready for use
016.015	Bacteria filter discs Ioniser (5 pcs.)
028.121	Grab rail ION Modul
026.118	Ionisation module external (round)
006.181	Power cable 4m
006.308	Chassis with black caps
006.188	Wheel for housing
006.348	Allen key for wheels 4mm

## 11 Disposal

medicap homecare GmbH complies with the criteria of the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU.

The product contains substances that may harm the environment if not disposed of in accordance with national legislation.

Therefore, the device or packaging can be returned to medicap homecare GmbH for disposal free of charge.

We ensure environmentally friendly disposal.

Do not dispose of used batteries with household waste!

Accessories that are not part of the product **MUST** be disposed of in accordance with the relevant disposal labelling.

**medicap homecare GmbH**

Hoherodskopfstr. 22

35327 Ulrichstein

Germany

Tel.: 06645/970-0

Fax: 06645/970-200