Aero life 2

Aero life 3

(DE) Gebrauchsanweisung	4
(GB) Instructions for Use	28
(FR) Mode d'emploi	50
(NL) Gebruiksaanwijzing	74
(IT) Istruzioni per l'uso	98

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Edition: Rev. P 2024-11
Date of issue: 2024-11-27

Table of contents

1 Introduction	30
1.1 Abbreviation.....	30
1.2 How to use the Instructions for Use.....	30
1.3 Significance of warnings and notes.....	30
1.4 Brief description.....	30
1.5 Intended use.....	31
1.6 Functional description.....	31
2. Safety and important information	32
2.1 Intended purpose.....	32
2.2 Indications and contraindications.....	32
2.2.1 Side effects and residual risks.....	33
2.3 Intended user and patient target groups.....	33
2.3.1 User responsibility.....	34
2.4 Important instructions for use and safety.....	34
2.4.1 General handling of the device.....	34
2.4.2 Handling recommended accessories.....	35
2.5 Warnings.....	35
2.6 Expected service life.....	37
2.7 Disclaimer of liability.....	37
2.8 Guarantee and warranty.....	37
2.9 Adresses.....	38
3 Operation	39
3.1 Preparation for operation.....	39
3.2 Switching ON.....	40
3.3 Operation.....	41
3.4 Switching OFF.....	41
4 Malfunctions and their elimination	42
5 Cleaning Process	43
5.1 General Notes.....	43
5.2 Cleaning Intervals.....	44
5.3 Cleaning and Disinfection.....	44
6 Maintenance	45
7 Pictorial Symbols	46
8 Technical Data and Specifications	47
9 Storage and Scope of delivery	48
10 Recommended Accessories, Disposables and Spare Parts	48
11 Disposal	49

1 Introduction

1.1 Abbreviation

MD	Medical Device
MDR	Medical Device Regulation (2017/745)
MEE	Medical electrical equipment
IPPB	Intermittend Positive Pressure Breathing

1.2 How to use the Instructions for Use

Importance of the Instruction for Use This Instructions for Use is part of the accompanying documents and are an essential part of the device. It contains all necessary instructions for operating the device in accordance with Regulation (EU) 2017/745, Annex I: General safety and performance requirements, Chapter III: Requirements regarding the information supplied with the device.

The Instruction for Use must be read thoroughly before the device is put into operation.

The device may only be used by persons who have been instructed in its use.

Illustrations The illustrations used in the documents may differ from the original.

Changes Changes to the documents are made as new editions or supplementary sheets. In general, the following applies: Subject to change.

Reproduction Reproduction, even in part, is only permitted with written approval.

1.3 Significance of warnings and notes

Warning

The signal word "Warning" alerts the user to a potentially hazardous situation which, if not avoided, could result in death or serious injury.

Caution

The signal word "Caution" draws the user's attention to the fact that in case of failure to observe this warning may result in injury to persons, serious injury or damage to the device.

Note

The signal word "Note" draws the user's attention to the fact that the following effects are to be expected in case of non-observance:

- Damage to the device may occur.
- Desired functions do not run or do not run correctly.

1.4 Brief description

The Aero life 2 / Aero life 3 excess pressure inhalator allows a patient with impaired or restricted breathing to breathe slowly and deeply with minimal effort.

The Aero life 2 / Aero life 3 helps the patient to breathe regularly by intermittent pressure build-up and enables simultaneous aerosol therapy.

Aero life 2

Aero life 3

- particularly suitable for home patients and practice

- operated by built-in, Low-noise compressor

- easily transportable

- for practice and clinic

- compressed air operated device, therefore be used without risk in potentially explosive atmospheres.

- universal installation (desk, wall, stand)

After surgical procedures, IPPB enables deep, slow breathing and better secretion clearance. The IPPB therapy, however, should only be carried out after thorough examination by a physician.

In case of pronounced pulmonary insufficiency, treatment must be carried out under the supervision of a physician.

1.5 Intended use

Precondition for the intended use of the device is the exact knowledge and observance of these Instructions for Use.

The device may only be used in accordance with the general regulations for the installation and operation of medical devices (MPBetreibV – Ordinance on Operators of Medical Devices).

Before using the device, the user must satisfy himself of its functional safety and proper condition.

Read all safety instructions and warnings carefully before using the device. Follow the safety instructions to avoid injuries and life-threatening situations.

When using the accessories, the Instruction for Use of the respective manufacturers and the corresponding warnings must be followed.

If additional devices are connected, the permissible leakage currents may be exceeded. Electromagnetic interference from outside does not pose a risk to the user.

1.6 Functional description

Intermittent positive pressure breathing (IPPB) is a periodical breathing of the lungs with air or oxygen assisted by slight pressure.

The Aero life 2 / Aero life 3 excess pressure inhalator allows a patient with impaired or restricted breathing to breathe slowly and deeply with minimal effort.

After an inspiration pulse from the patient, the Aero life 2 / Aero life 3 pumps an air-oxygen mixture into the patient's lungs at the previously set pressure. When the desired pressure is reached in the lungs, the device automatically switches to expiration - the patient can breathe out freely. Under normal conditions, the Aero life 2 / Aero life 3 automatically follows the patient's breathing rate.

The Aero life 2 / Aero life 3 helps the patient to breathe regularly by intermittent pressure build-up and enables simultaneous aerosol therapy.

To avoid damage to the sensitive components, the proper handling of the device and its accessories must be observed. For this purpose, all instructions and warnings, as well as the exact operation or preparation of the components and their cleaning and maintenance must be observed.

2. Safety and important information

2.1 Intended purpose

The Aero life 2 / Aero life 3 IPPB- device (Intermittent Positive Pressure Breathing) enables a patient with **an** impaired or restricted breathing to breathe slowly and deeply with minimal effort.

The operating environment must be complied in accordance with the environmental conditions specified in chapter 8 Technical data and Specifications. The warnings in chapter 2.5 must be observed.

Risks were evaluated in the risk management.

Caution

The physician's instructions must be followed strictly.

Caution

Any physical disorder must be reported to the attending physician immediately.

Caution

The effectiveness of the therapy must be checked regularly.

Warning

Improper handling and no intended use can cause danger and damage. Therefore, we ask you to read through these Instruction for Use and follow them strictly. Always keep them within easy reach.

2.2 Indications and contraindications

Indications

Application should only be carried out in accordance with the physician's instructions. The duration and frequency of use should depend on the disease.

The recognised treatment concept of IPPB therapy is mainly used in the treatment of chronic obstructive pulmonary diseases such as

- COPD (chronic obstructive pulmonary disease)
- Respiratory
- Pulmonary fibrosis, Pulmonary emphysema
- Pneumoconiosis
- Bronchitis, bronchial asthma and other infections or bronchospastic changes of the lungs

IPPB therapy is intended to help patients with impaired or restricted breathing to achieve regular deep ventilation.

The patient group consists of patients suffering from respiratory diseases. Diseases of the respiratory tract are often chronic and therefore often need to be treated several times a day.

Contraindications

Contraindications of the accessories can be found in the accompanying documents (e.g. Instructions for Use) of the accessories.

IPPB therapy is dependent on the following conditions, among others:

- Age of the patient
 - General condition of the patient (e.g. fever)
 - Smoking or coffee consumption
- (This list is not complete, but merely exemplary).

It may only be used on high-risk patients with particular caution and under the special supervision of the attending physician.

High-risk patients are persons with:

- Advanced age
- Cerebral spasms
- Fever
- Severe lung diseases

Caution

Any physical disorder must be reported to the attending physician immediately.

Caution

The IPPB- device (Intermittent Positive Pressure Breathing) is not intended for life-supporting or life-sustaining use.

2.2.1 Side effects and residual risks

The following side effects may occur as a result of interactions if the device is used improperly:

- Overdosage of drugs/ medicine.

The side effects and contraindications of the medications used must be observed.

2.3 Intended user and patient target groups

Patient target groups

The intended target group are patients with the above-mentioned indications.

The physician's instructions must be followed strictly.

Permissible age, gender, duration and frequency of use is to be determined by the attending physician.

The application is not intended for patients with contraindications. These belong to the patient risk groups.

Caution

The device is not intended for life-supporting or life-sustaining use.

Warning

Any physical disorder must be reported to the attending physician immediately.

User group

The device may only be set up, operated and used by persons who have the required (necessary) training or knowledge and experience and who have been successfully instructed in its use

The IPPB- devices (Intermittent Positive Pressure Breathing) are used in domestic, commercial and clinical settings.

The users are:

- Physicians
- Alternative practitioner
- Auxiliary / carer
- Private persons

The Aero Life 2 IPPB- device (Intermittent Positive Pressure Breathing) is particularly suitable for home patients and surgeries.

The Aero Life 3 IPPB- device (Intermittent Positive Pressure Breathing) is a compressed air-operated device and therefore suitable for surgeries and clinics.

Patients under the age of 18 or patients with impaired mental or physical condition will only be treated under the supervision of another qualified person.

Patients under 6 years of age may only be treated under the supervision of another qualified person.

Application /therapy in pregnant or breastfeeding women should only be carried out in accordance with the physician's precise instructions

2.3.1 User responsibility

The device may only be operated /used under the operating conditions specified by the manufacturer.

The user may only use the device in the language in which the Instruction for Use is provided.

Note the following when entering parameters:

- The parameters entered must be verified by the user, i.e. the user must check the correctness of the values entered.

The user and/or patient must report all serious incidents related to this device to the manufacturer and/or the competent authority of the Member State in which the user and/or patient is established.

The user does not need any specific qualification, but must follow the requirements for the intended use of the device, depending on the manufacturer's specifications, in these Instruction for Use

The device may only be used by persons who have been instructed in its use.

2.4 Important instructions for use and safety

2.4.1 General handling of the device

Note

Precondition for the intended use of the device is the exact knowledge and observance of this Instruction for Use.

Note

The device must be installed in a well-ventilated, dry, cool and dust-free room.

Note

The information on the identification plate and the instructions on the labels must be observed.

Note

If the device has been stored at temperatures below 10 °C, it should be allowed to adjust to room temperature for approx. 60 minutes, otherwise malfunctions may occur.

Note

Keep the device free from oil and grease.

Note

Only original parts should be used. Do not use worn or damaged parts, such as the mains cable.

Note

The device must only be opened by an authorized service technician. No changes or modifications may be made to the device.

Note

Do not connect the device to the supply voltage via a multiple socket. The mains cable must not exceed a maximum length of 4m.

Caution

The device must be protected from moisture and humidity. DO NOT touch the device when it is wet.

Caution

There is no danger to the user from external electromagnetic interference.

Note

Check ALL external components and the packaging for damage. In the case of damage or if the product does not function properly, contact a service technician or the manufacturer so that appropriate repair measures can be carried out.

Note

If possible, the device should not be installed near radiators, fan heaters, etc.

Caution

The air supply must not be impaired, as cooling air is required for the compressor. The minimum distance from cabinets etc. should be 60 cm.

Caution

Maintenance and servicing must only be carried out by an authorized service technician.

Caution

The device must be checked annually by a person authorized by medicap homecare GmbH.

Caution

When not in use for a longer period of time, unplug the power supply from the socket.

2.4.2 Handling recommended accessories

Caution

The specified accessories, disposables and spare parts have been validated by medicap homecare GmbH for use with the devices.
When using accessories, disposables and spare parts from other manufacturers, accuracy can no longer be guaranteed.

Note

The accessories, which are declared as "single use" article, are not suitable for repeated use.

Caution

Reusing the single-use disposables, such as the Tubing System IPPB System ML, can endanger both the patient and the operating personnel.

Caution

During medication nebulisation, only the medication specified by the physician may be used.

2.5 Warnings

Warning

The device may only be operated by the user under the operating conditions specified by the manufacturer, see chapter 8 Technical data.

Warning

Compliance with all instruction for use and safety instructions from Chapter 2 is mandatory.

Warning

Smoking, open light and fire are prohibited during use.

Warning

The device must not be used if:

- the mains plug or cable is damaged
- functions run incorrectly
- the device has been dropped
- the device is damaged
- liquid has penetrated the device

In these cases, the device must be checked by an authorized service technician and repaired if necessary.

Warning

The device must only be opened by an authorised service technician.

Warning

Multiple sockets or extension cables must not be used.

Warning

The device should not be placed near radiators, fan heaters, etc.

Warning

Portable and mobile HF communication equipment can affect medical electrical devices. This can lead to malfunctions in the devices.

For this reason, radiation sources (devices that emit electromagnetic waves) must be kept at a certain minimum distance from operating medical electrical devices. See table below:

Radiation sources

Minimum distance from a medical device including all connection cables

Mobile phone	3.3 m
DECT telephone (cordless telephone)	1.2 m
Bluetooth devices (laptops, mobile phones)	0.7 m
Remote controls	1.2 m
WLAN devices (e.g.: laptops, repeater, access point, print server)	2.3 m

Warning

DO NOT use the device in an environment where there is contamination, smoke, gases, flammable anesthetics, cleaning agents or chemical vapors.

Warning

The use of accessories and cables other than those listed in the Instruction for Use can lead to increased interference emission, reduced interference immunity of the device or impairment of the application.

Warning

Close supervision is required when this device is used in the vicinity of children or disabled persons.

Warning

Patients who cannot see or hear alarm signals or cannot communicate their discomfort must be supervised.

Warning

A change in altitude above sea level can affect the performance of the device. Consult your physician or contact the manufacturer.

Warning

Stow and lay electrical cables properly to avoid potential tripping hazards.

Warning

Risk of infection if the device is used again!
If the device is used by several patients, infections can be transmitted to the next patient.

When the device is used again, it must be cleaned and disinfected, see chapter 5.

Warning

The ME equipment has no defibrillation protection.

Warning

Modification of the ME equipment is not permitted.

Warning

The device may not be modified without the manufacturer's permission.

Warning

If the ME equipment is modified, suitable inspections and tests must be carried out to ensure continued safe use, suitable inspections and tests must be carried out.

Warning

If the mains supply is interrupted, the ME equipment has no power supply.

2.6 Expected service life

If maintenance is carried out to the prescribed extent and at the specified intervals, the device will operate safely during the intervening period.

With each maintenance and Technical Safety Control, the “expected service life” according to IEC 60601-1 will therefore be prolonged until the next prescribed Technical Safety Check.

2.7 Disclaimer of liability

The device has been approved for use with certain accessories and spare parts (see Chapter 10).

The manufacturer accepts no responsibility or liability for personal injury or other damage and excludes any warranty for damage to the device resulting from the use of non-approved or unsuitable accessories or spare parts

2.8 Guarantee and warranty

medicap homecare GmbH provides a two-year guarantee from the date of delivery for defects caused by material or manufacturing faults.

Defects covered by the warranty will be rectified in accordance with the warranty conditions.

In addition, medicap homecare GmbH does not provide a guarantee if the operator jeopardizes the functions of the device by failing to observe these Instructions for Use, by using the device for purposes other than those for which it is intended or by tampering with it.

medicap homecare GmbH does not automatically consider the owner of the device to be authorised to carry out service work.

Note

The guarantee can only be claimed in conjunction with the proof of purchase.

2.9 Adresses

Manufacturer:

medicap homecare GmbH
Hoherodskopfstr. 22
35327 Ulrichstein
Germany

Tel.: +49 6645/970-0
Fax: +49 6645/970-200


3 Operation

3.1 Preparation for operation

The Aero life 2 / Aero life 3 must be set up in a well-ventilated, dry, cool and dust-free room. The air supply must not be impaired, as cooling air is needed for the compressor.

- With the Aero life 2, the mains plug must be plugged into a socket.
- With the Aero life 3, the compressed air plug must be connected to a compressed air system.

The Tubing System IPPB System ML must be connected. The tubings have three different diameters and only fit on the designated connections of the same size on the device.

Caution

The settings are to be carried out according to the physician instructions.

Delivery setting:

	<u>delivery setting Aero life 2</u>	<u>delivery setting Aero life 3</u>
Ventilation pressure	15- 20	15
Flow rate	30 l/ min	30 l/ min
Starting threshold	-1,5 mbar	-1,5 mbar
Nebuliser output	maximum setting	maximum setting
Exhalation resistance	0	0


3.2 Switching ON

If the device has been stored at temperatures below 10 °C, the device must adjust to room temperature, otherwise malfunctions may occur.

→ Aero life 2: Press the ON/OFF button, the compressor switches on and the device is ready for operation.

→ Aero life 3: Press the toggle switch and the device is ready for operation.

3.3 Operation

The medication nebuliser must be used horizontally.

Caution

During medication nebulisation, only the medication specified by the physician may be used.

Place the mouth piece in the mouth. The suction on the mouthpiece triggers the inhalation. Breathing is now supported by the device. When the set positive pressure is reached, the device switches to exhalation.

Caution

Ensure that no air escapes through the nose (use nose clamps if necessary).

3.4 Switching OFF

→ Aero life 2: Press the ON/OFF button.

→ Aero life 3: Press the toggle switch.

If the device is not to be used for a longer period of time, the plug must be pulled out of the socket (Aero Life 2) or disconnected from the reducing valve or the compressed air system (Aero life 3).

The Tubing System IPPB System ML is intended for single use only and must be disposed of after use.

4 Malfunctions and their elimination

Failure symptom	Cause of error	Remedy
Compressor does not run	Mains plug not inserted into the socket	Insert the mains plug into the socket
	Mains fuse defective	Change (replace) mains fuse 0.8 AT
Device does not build up pressure	Incorrect setting	Carry out the settings specified by the physician or restore the delivery setting (see chapter 3.1 Preparing for operation).
Start threshold not reached (device does not start up on inhalation)	Tubing System IPPB System ML not connected correctly	Connect Tubing System IPPB System ML correctly
	Expiration valve not installed correctly (Tubing System IPPB System ML)	Connect the tube correctly to the elbow of the expiration valve
	Controller start threshold not correctly set	Carry out the settings specified by the physician or restore the delivery setting (see chapter 3.1 Preparing for operation).
No nebulisation	Nebuliser output set too low	Carry out the settings specified by the physician or restore the delivery setting (see chapter 3.1 Preparing for operation).
	Tubing System IPPB System ML not connected correctly	Connect Tubing System IPPB System ML correctly
	Venturi nozzle in the medication container stuck with medication (Tubing System IPPB System ML)	Replace the Tubing System IPPB System ML or clean the Venturi nozzle.
Exhalation more difficult	Medication container empty	Fill medication container
	Exhalation resistance set incorrectly (too high)	Carry out the settings specified by the physician or restore the delivery setting (see chapter 3.1 Preparing for operation).

5 Cleaning Process

Equipment cleaning is a process that a used product undergoes so that it can be safely reused. This procedure includes cleaning and disinfection.

In addition, a test and restoration of the technical and functional safety of the product used must be performed carried out.

Warning

Risk of infection if the device is used again! If the device is used by several patients, infections can be transmitted to the next patient. If the device is used again, it must be cleaned and disinfected.

5.1 General Notes

Note

Wear suitable protective equipment during cleaning / disinfection.

Caution

The cleaning and disinfecting agent must be used in accordance with the manufacturer's instructions, see the manufacturer's Instructions for Use.

Caution

Do not allow any liquid to enter the interior of the device.

Caution

DO NOT allow detergent to drip into the air inlet or outlet ports.

Caution

After cleaning, disinfection and a functional check, the device is suitable for reuse with further patients.

Warning

Risk of injury due to electric shock! Penetrating liquids from cleaning agents and disinfectants can cause a short circuit, injure the user and damage the device.

Caution

Ensure spatial separation of equipment that is still contaminated and equipment that has already been disinfected (such as work surfaces and storage rooms, etc.).

Caution

Switch off the device and disconnect it from the power supply before starting cleaning.

Caution

Do not open the device! Repairs may only be carried out by qualified technicians.

Caution

Do not pour liquid over the device and components or immerse them in liquids.

Caution

DO NOT clean the housing with concentrated chlorine-based products (ethylene chloride), oil-based products (Pine-Sol®, Lestoil®), or other aggressive chemical cleaning agents. Use commercially available cleaning agents to clean the outer housing, but do not use abrasive, scratching or corrosive agents.

Caution

No lubricants other than those recommended by medicap homecare GmbH may be used.

5.2 Cleaning Intervals

Intervals	Action
daily	<ul style="list-style-type: none"> - Clean / disinfect flow block / patient connection - Accessories which are declare as "single use" are not for repeated use
weekly	<ul style="list-style-type: none"> - Clean / disinfect housing - Clean / disinfect operating elements
monthly	
annual	- the device must be checked annually by personnel authorized by medicap homecare GmbH as part of the Technical Safety Control.
when needed	<ul style="list-style-type: none"> - Clean / disinfect housing - Clean / disinfect operating elements - Clean / disinfect flow block / patient connection - Clean / disinfect mains cable - Replace / Change approved disposables
For patient change	<ul style="list-style-type: none"> - Clean / disinfect housing - Clean / disinfect flow block / patient connection - Clean / disinfect all operating elements - replace / exchange approved disposables

5.3 Cleaning and Disinfection

Component / part	Cleaning (manual)	Disinfection (manual)
Housing Operating elements Flow block / Patient connection	Wipe with a damp cloth; Use water or commercially available cleaning agents, but do not use abrasive, scratching or corrosive agents, allow to dry.	Wipe disinfection (see recommended disinfectants below), allow to dry.
Power (Mains) cable	Wipe with a damp cloth; Use water or commercially available cleaning agents, allow to dry	Wipe disinfection (see recommended disinfectants below), allow to dry.
Accessories such as the Tubing System IPPB System ML	Single-use items are to be disposed of, no cleaning	Single-use items are to be disposed of, no disinfection

Cleaning of the devices is always done manually. The cleaning agent is to be applied on a cloth (cloths with low particle release are to be used).

Brushes, spray guns or similar are not recommended.

Recommendation for cleaning agents and disinfectants

The following disinfectants have been validated by medicap homecare GmbH and can be used:

- Biguanid surface N (Dr. Schuhmacher)
- Descogen Liquid
- Korsolex (Bode Chemie)
- mikrozyd AF wipes (Schülke)

Any other disinfectant used must be listed in the VAH list and be compatible with the materials to be treated.

Drying is carried out in fresh air without certain drying media (such as hair dryers). Wiping with cloths (for drying) is not recommended.

6 Maintenance

Caution

All maintenance intervals are guide values. In individual cases, these may be shorter due to changes in environmental conditions with regard to temperature, humidity, dust, etc.

Caution

Maintenance may only be carried out by an authorized service technician.

Caution

Only the original filter may be used.

Caution

Circuit diagrams, component lists, descriptions or other information can be provided on request.

Caution

During maintenance, only monitored test equipment and measuring devices suitable for the intended purpose may be used.

Regular maintenance is required once a year, as wearing parts are built into the device.

With regular use, the check valve flow block, the double silencer and the built-in bacteria filter must be serviced.

The device must not be maintained or serviced during use on the patient!

Caution

Works on the opened device may only be performed by an authorized service technician.

7 Pictorial Symbols

Attention!

Smoking is forbidden

Follow the Instructions for Use

protection class II

Single-use; Do not reuse

Serial number

The device must not be disposed of with household waste

CE- Mark 0044; Notified body: TÜV NORD CERT GmbH 0044

Degree of protection against electric shock: Type BF

Alternating current

Dangerous electric-voltages

Manufacturer

Date of manufacture

Medical device

UDI number - Unique device identifier

8 Technical Data and Specifications

Dimensions:	23 x 15 x 27 cm
Weight:	approx. 5 kg
Ambient temperature Operating:	10 to 40 °C
Operating voltage:	230 V+/-10% 50 Hz
Power consumption:	approx. 65 W
Energy consumption:	approx. 0,065 kWh
Pressure adjustment:	0 – 35 mbar
Flow rate:	20 – 60 l/min
Starting threshold:	-0,5 bis – 5mbar
Nebuliser capacity:	infinitely variable
Exhalation resistance:	infinitely variable
Control:	assistiert
Protection class:	II
Degree of protection against electric shock:	Type BF
MDD 93/42 EEC classification:	Ila
Device corresponds:	MDD 93/42/EEC
Relative humidity:	non-condensing 10 to 95 %
Storage temperature:	-10 to 50°C
Degree of protection against penetration of liquids:	IP 22
Notified body:	 0044, TÜV Nord Cert

9 Storage and Scope of delivery

The device is delivered safely packed in the box. Transport and storage take place in the original box. Store the device, before and after use, under the prescribed environmental conditions. If the device has been used, the device should be cleaned and disinfected before storage.

Description	Article number
Scope of delivery Aero life 2, consisting of: <ul style="list-style-type: none"> - Aero life 2 - Tubing System IPPB System ML - Holder IPPB-System ML (Aero life 2) - Instruction for Use - Bag Aero life 2 black 	002.000
Scope of delivery Aero life 3, consisting of: <ul style="list-style-type: none"> - Aero life 3 - Tubing System IPPB System ML - Holder IPPB-System ML - Instruction for Use - PVC Fabric tubing 6x3mm neutral 	003.000

10 Recommended Accessories, Disposables and Spare Parts

The accessories, disposables or spare parts listed below have been validated by medicap homecare GmbH for use with the Aero Life 2/3 medical device and are recommended.

Caution

The accessories such as the Tubing System IPPB System ML are "single use" items. They are NOT suitable for repeated use by a single patient as well as by multiple patients in the physician's office/clinic.

Caution

The use of accessories, spare parts or disposables that have NOT been approved may result in impaired performance.

Caution

Only the original filter may be used.

Article or spare part number	Description
002.010	Tubing System IPPB System ML
002.078	Holder for Tubing System IPPB System ML
002.198	Instruction for Use Aero life 2 / Aero life 3
002.017	Bag Aero life 2 black
002.043	Trolley for Aero life 2 / Aero life 3)
002.146	Mains cable for Aero life 2 with Euro central plug white

11 Disposal

medicap homecare GmbH complies with the criteria of the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU.

The product contains substances that may harm the environment if not disposed of in accordance with national legislation.

Therefore, the device or packaging can be returned to medicap homecare GmbH for disposal free of charge. We will ensure that it is disposed of in an environmentally friendly manner.

The disposal of accessories which are not part of the product **MUST** be carried out in accordance with the respective disposal label.

medicap homecare GmbH
Hoherodskopfstr. 22
35327 Ulrichstein
Germany
Tel.: + 49 6645/970-0
Fax: +49 6645/970-200

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