medicap
homeCare GmbH

Précise 8008

(DE) Gebrauchsanweisung	3
(GB) Instructions for Use	52
(ES) Instrucciones de uso	102
(IT) Istruzioni per l'uso.....	152

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

Display Patientendaten ändern:

The screenshot shows a user interface for changing patient data. At the top, there are four circular icons: a person with a green cross, a person with a pencil, a person with a chart, and a blue square with a white arrow pointing right. The second icon (person with a pencil) is highlighted with a red box. Below these are four input fields: a placeholder for 'Vorname' with a person icon, a placeholder for 'Nachname' with a person icon, a date input field showing '12 . 09 . 1980' with a calendar icon, and a gender selection field with a male icon. A large circular icon with a person and a pencil is centered at the bottom.

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Display Patientendatenbank:

The screenshot shows a user interface for displaying patient data. At the top, there are four circular icons: a person with a green cross, a person with a pencil, a person with a chart, and a blue square with a white arrow pointing right. The second icon (person with a pencil) is highlighted with a red box. Below these are two input fields: one for 'Vorname' with a person icon and one for 'Nachname' with a person icon. A table below lists patient data: Hans (Vorname), Mustermann (Nachname), 12.09.1980 (Geburtsdatum), and a male icon (Geschlecht). A large circular icon with a person and a pencil is centered at the bottom.

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
Display Systemeinstellungen:

The screenshot shows a vertical list of icons on the left and a detailed configuration table on the right. The icons are: information (i), trash can, printer with gear, volume, and a blue square with a white arrow pointing right. The table has two columns: 'Aktive Eventliste' and 'E1; E2; E3; E4; E5'. It lists software versions, battery status, controller info, sensor versions, and a language selection dropdown. The language dropdown shows flags for German, English, French, Czech, Polish, Spanish, Portuguese, and Greek, with 'Deutsch' selected. An arrow points from the text 'Sprachauswahl' to this dropdown.

Aktive Eventliste	E1; E2; E3; E4; E5
Software	v3.2 (11.05.2022)
DB-Texte	v2.0.0 (28.03.2024)
DB-Messdaten	v2.0 (07.02.2020)
Akku-Stand	97 %
Akku-Spannung	15.5 V
Controller	Feb 21 2020
Sensor1	v0.13
Sensor2	v0.13
Sensor3	v0.13
Sensor4	v0.13
Sensor5	v--
Sensor6	
Sensor7	
Sensor8	
Geräte-ID	0.0100
Geräte-IP	
Standarddrucker	0-series
Sprache	

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

The figure consists of two vertically stacked screenshots of a medical software interface. On the left is a vertical toolbar with icons for information, trash, printer, volume, and a clipboard. The main area shows patient data in a table.

Screenshot 1 (Top): A table with columns Vorname, Nachname, and Geburtsdatum. The row contains Hans, Mustermann, and 12.09.1980. To the right of the table are three icons: a trash can, a megaphone, and a circular arrow. An arrow points from the text "Patient auswählen" to the trash can icon.

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Screenshot 2 (Bottom): The same table structure, but the "Geburtsdatum" column has a checked checkbox in its header. Below the table is another table with columns Datum and Uhrzeit, containing rows for 02.04.2024, 14:38:59; 28.03.2024, 14:07:20; and 28.03.2024, 12:34:52. To the right of this second table are four icons: a trash can, a megaphone, a circular arrow, and a checkmark. Arrows point from the following labels to these icons:

- "Patient ist ausgewählt" to the checkmark icon above the first table.
- "Patient löschen" to the trash can icon next to the first table.
- "Auswahl aufheben" to the circular arrow icon next to the first table.
- "Messung auswählen" to the checkmark icon above the second table.
- "Messung löschen" to the trash can icon next to the second table.
- "Auswahl aufheben" to the circular arrow icon next to the second table.

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**Edition: Rev. C 2024-10
Date of issue: 2024-10-01**

Table of contents

1 Introduction	55
1.1 Abbreviation	55
1.2 How to use the Instructions for Use	55
1.3 Significance of warnings and notes	55
1.4 Brief description.....	55
1.5 Intended use.....	56
1.6 Functional description	56
2 Safety and important information	57
2.1 Intended purpose.....	57
2.2 Indications and contraindications.....	57
2.2.1 Side effects and residual risks	58
2.3 Intended user and patient target groups.....	58
2.3.1 User responsibility	59
2.4 Important instructions for use and safety	59
2.4.1 General handling of the device	59
2.4.2 Handling recommended accessories.....	60
2.5 Warnings	60
2.6 Expected service life	62
2.7 Disclaimer of liability.....	63
2.8 Guarantee / Warranty	63
2.9 Adresses.....	63
3 Operation	64
3.1 Preparation	64
3.1.1 Connecting the power supply unit.....	64
3.1.2 Assignment of the sensors	64
3.1.3 Connecting sensors	65
3.1.4 Fixation & use of the supporting frame with a retaining arm.....	66
3.2 Operation	67
3.2.1 General	67
3.2.2 Switching on the device	67
3.2.3 Fixing the sensor and preparing the measurement site	68
3.2.4 General operation - Patient menu and patient database	69
3.2.5 Measurement	72
3.2.6 Sensor removing	76
3.2.7 Evaluate measurement data.....	77
3.2.8 System menu.....	79
3.2.8.1 System settings.....	79
3.2.8.2 Delete data	81
3.2.9 Date & time.....	82
3.2.10 Switching off the device.....	83
3.3 Print functions	84
3.3.1 USB interface.....	84
3.3.2 Define settings	85
3.3.3 Print measurement	86
3.3.4 Data export	88
3.4 Sensor menu.....	90

3.4.1 Set sensor temperature.....	90
3.4.2 Switching the sensor on/off	91
4 Malfunctions and their elimination	92
5 Cleaning Process	92
5.1 General Notes.....	92
5.2 Cleaning Intervals	93
5.3 Cleaning and Disinfection.....	94
6 Maintenance.....	95
7 Pictorial Symbols	96
7.1 Definition of the symbols on the device or accessories	96
7.2 Software symbols.....	97
8 Technical Data and Specifications.....	99
9 Storage and Scope of delivery	100
10 Recommended Accessories, Disposables and Spare Parts	100
11 Disposal	101

1 Introduction

1.1 Abbreviation

MD	Medical Device
MDR	Medical Device Regulation (2017/745)
ME	Medical electrical equipment

1.2 How to use the Instructions for Use

Importance of the Instruction for Use This Instructions for Use is part of the accompanying documents and therefore an essential part of the device. It contains all necessary instructions for operating the device in accordance with Regulation (EU) 2017/745, Annex I: General safety and performance requirements, Chapter III: Requirements regarding the information supplied with the device.

The Instruction for Use must be read thoroughly before the device is put into operation.

The device may only be used by persons who have demonstrably been instructed in its proper use.

Illustrations The illustrations used in the documents may differ from the original.

Changes Changes to the documents are made as new editions or supplementary sheets. In general, the following applies: Subject to change.

Reproduction Reproduction, even in part, is only permitted with written approval.

1.3 Significance of warnings and notes

Warning

The signal word "Warning" alerts the user to a potentially hazardous situation which, if not avoided, could result in death or serious injury.

Caution

The signal word "Caution" draws the user's attention to the fact that in case of failure to observe this warning may result in personal injury, serious injury or damage to the device.

Note

The signal word "Note" draws the user's attention to the fact that the following effects are to be expected in case of non-observance:

- Damage to the device may occur.
- Desired functions do not run or do not run correctly.

1.4 Brief description

The Précise 8008 transcutaneous oxygen measuring device is used to measure the partial pressure of oxygen, also known as tcpO₂, on the surface of the skin. The measured values can be displayed graphically, saved and printed out.

The sensor parameters, the measurement process and the measurement results can be saved in the internal database. This means that every measurement can be reliably documented and traced at any time.

Moreover, visualization software (oder visual display software) can also be used, which offers the user a further option for graphically displaying the recorded measurement data in addition to the Précise 8008 display. In addition to the extended display, there are other functions that are explained in detail in this document.

This results in two versions of the Précise 8008:

Version 1:

The computer is connected to the interface on the side of the Précise 8008 using a USB cable. The data can be displayed and saved in real time with the visualization software via this connection. Furthermore, a detailed tcpO2 measurement report can be saved as a PDF document or printed out. A pure "raw data export" is also possible.

Version 2:

The Précise 8008 is used autonomously as a "standalone device". Using a USB stick, the measurement data can be recorded, called up and displayed in a "raw data form". This means that every measurement can be documented and processed externally on the PC at any time.

Caution

With the interface box, the device is also suitable for measurement in a pressure chamber. The device must NOT be in a pressure chamber, only the interface box.

Warning

The device must not be operated in a pressure chamber. If a measurement is carried out in a pressure chamber, only the sensor and the accessories, such as the interface box, may be in the pressure chamber.

1.5 Intended use

Precondition for the intended use of the device is the exact knowledge and observance of these Instructions for Use.

The device may only be used in accordance with the general regulations for the installation and operation of medical devices (MPBetreibV – Ordinance on Operators of Medical Devices).

Before using the device, the user must satisfy himself of its functional safety and proper condition.

Read all safety instructions and warnings carefully before using the device. Follow the safety instructions to avoid injuries and life-threatening situations.

When using the accessories, the Instruction for Use of the respective manufacturers and the corresponding warnings must be followed.

If additional devices are connected, the permissible leakage currents may be exceeded. Electromagnetic interference from outside does not pose a risk to the user.

1.6 Functional description

The tcpO2 measurement is a non-invasive method for determining the partial pressure of oxygen on the skin surface. During the procedure, the skin is warmed to 40° - 44°C.

The oxygen measurement is based on the luminescence lifetime measurement. This is dependent on the oxygen partial pressure. Measurement preparations (e.g. changing the electrolyte, covering the sensor, etc.) are no longer necessary, as is known from the Clark sensor.

The sensors are suitable for both the determination and long-term measurement of oxygen partial pressure.

To avoid damage to the sensitive components, the proper handling of the device, sensors and its accessories must be observed. For this purpose, all instructions and warnings, as well as the exact operation or preparation of the components and their cleaning and maintenance must be observed.

2 Safety and important information

2.1 Intended purpose

The Précise 8008 transcutaneous oxygen measuring device is used to measure the partial pressure of oxygen, also known as tcpO₂, non-invasively on the surface of the skin. The measurement values can be displayed graphically, saved and printed out.

The operating environment must be complied in accordance with the environmental conditions specified in chapter 8 Technical data and Specifications. The warnings in chapter 2.5 must be observed.

Risks were evaluated in the risk management.

Caution

The physician's instructions must be followed strictly.

Caution

Any physical disorder must be reported to the attending physician immediately.

Warning

Improper handling and no intended use can cause danger and damage. Therefore, we ask you to read through these Instruction for Use and follow them strictly. Always keep them within easy reach.

2.2 Indications and contraindications

Indications

Application should only be carried out in accordance with the physician's instructions. The duration and frequency of use should depend on the disease.

Current areas of application include the assessment of peripheral arterial occlusive disease, chronic venous insufficiency and diabetes mellitus as well as collagenoses and inflammatory skin diseases.

In addition to phlebological / angiological microcirculation diagnostics, clinical areas of application include anesthesia and intensive care medicine as well as respiratory physiology and sports medicine.

Contraindications

The use of this device is unsuitable for various skin hypersensitivities.

Contraindications of the accessories can be found in the accompanying documents (e.g. Instructions for Use) of the accessories.

The measured transcutaneous oxygen partial pressure is dependent on the following conditions, among others:

- Temperature selection
- Measurement site selection
- Age of the patient
- General condition of the patient (e.g. fever)
- Smoking or coffee consumption
- Acclimating patients to the outside temperature
(This list is not complete, but merely exemplary).

It may only be used on high-risk patients with particular caution and under the special supervision of the attending physician.

High-risk patients are persons with:

- Advanced age
- Cerebral spasms
- Viral infectious diseases
- Epilepsy
- Somatic diseases
- Simultaneous ACTH or glucocorticoid treatment
- High carbon dioxide concentration in arterial (oxygen-rich) blood
- Poisoning by substances that reduce respiratory activity
- Fever

2.2.1 Side effects and residual risks

If the contraindications are observed, no side effects are to be expected when the measuring device is used as intended.

The following side effects may occur as a result of interactions if the measuring device is used improperly:

- Reddening of the skin

2.3 Intended user and patient target groups

Patient target group

The intended target group are patients with the above-mentioned indications.

The device is only intended for patients for whom the partial pressure of oxygen, also known as tcpO_2 , is to be measured at the skin surface.

The physician's instructions must be followed strictly.

Permissible age, gender, duration and frequency of use is to be determined by the attending physician.

The application is not intended for patients with contraindications. These belong to the patient risk groups.

Transcutaneous oxygen partial pressure measurement can also be carried out in patient risk groups (e.g. elderly patients) ONLY under the special care and supervision of the attending physician.

The resulting biological and toxicological risks must be checked and evaluated in advance by the attending physician.

Caution

The device is not intended for life-supporting or life-sustaining use.

Warning

Any physical disorder must be reported to the attending physician immediately.

User group

Clinical use

The device may only be set up, operated and used by persons who have the required (necessary) training or knowledge and experience and who have been successfully instructed in its use.

The users are:

- Physicians, Alternative practitioner, Auxiliary / carer

Patients under the age of 18 or patients with impaired mental or physical condition will only be treated under the supervision of another qualified person.

Patients under 6 years of age may only be treated under the supervision of another qualified person.

Application / measuring in pregnant or breastfeeding women should only be carried out in accordance with the physician's precise instructions

2.3.1 User responsibility

The device may only be operated /used under the operating conditions specified by the manufacturer.

The user may only use the device in the language in which the Instruction for Use is provided.

Note the following when entering parameters:

- The parameters entered must be verified by the user, i.e. the user must check the correctness of the values entered.

The user and/or patient must report all serious incidents related to this device to the manufacturer and/or the competent authority of the Member State in which the user and/or patient is established.

The user does not need any specific qualification but must follow the requirements for the intended use of the device, depending on the manufacturer's specifications, in these Instruction for Use.

The device may only be used by persons who have been instructed in its use.

2.4 Important instructions for use and safety

2.4.1 General handling of the device

Note

Precondition for the intended use of the device is the exact knowledge and observance of this Instruction for Use.

Note

The device must be installed in a well-ventilated, dry, cool and dust-free room.

Note

The information on the identification plate and the instructions on the labels must be observed.

Note

If the device has been stored at temperatures below 10 °C, it should be allowed to adjust to room temperature for approx. 60 minutes, otherwise malfunctions may occur.

Note

Keep the device free from oil and grease.

Note

Only original parts should be used. Do not use worn or damaged parts, such as the mains cable.

Note

The device must only be opened by an authorized service technician. No changes or modifications may be made to the device.

Note

Do not connect the device to the supply voltage via a multiple socket. The mains cable must not exceed a maximum length of 4 m.

Caution

The device must be protected from moisture and humidity. DO NOT touch the device when it is wet.

Caution

There is no danger to the user from external electromagnetic interference.

Note

Check ALL external components and the packaging for damage. In the case of damage or if the product does not function properly, contact a service technician or the manufacturer so that appropriate repair measures can be carried out.

Note

If possible, the device should not be installed near radiators, fan heaters, etc.

Caution

Maintenance and servicing must only be carried out by an authorized service technician.

Caution

The device must be checked annually by a person authorized by medicap homecare GmbH.

Caution

When not in use for a longer period of time, unplug the power supply from the socket.

2.4.2 Handling recommended accessories

Caution

The specified accessories, disposables and spare parts have been validated by medicap homecare GmbH for use with the devices.

When using accessories, disposables and spare parts from other manufacturers, accuracy can no longer be guaranteed.

Note

The accessories, such as the adhesive rings, are "single use" items. This means that they are NOT suitable for repeated use by a single patient or multiple patients.

Caution

Reusing the single-use disposables, such as the adhesive rings, can endanger both the patient and the operating personnel.

2.5 Warnings

Warning

The device may only be operated by the user under the operating conditions specified by the manufacturer, see chapter 8 Technical data.

Warning

Compliance with all instruction for use and safety instructions from Chapter 2 is mandatory.

Warning

Smoking, open light and fire are prohibited during use.

Warning

The device must not be used if:

- the mains plug or cable is damaged
- functions run incorrectly
- the device has been dropped
- the device is damaged
- liquid has penetrated the device

In these cases, the device must be checked by an authorized service technician and repaired if necessary.

Warning

The device must only be opened by an authorized service technician.

Warning

Multiple sockets or extension cables must not be used.

Warning

The device should not be placed near radiators, fan heaters, etc.

Warning

Portable and mobile HF communication equipment can affect medical electrical devices. This can lead to malfunctions in the devices.

For this reason, radiation sources (devices that emit electromagnetic waves) must be kept at a certain minimum distance from operating medical electrical devices. See table below:

Radiation sources

**Minimum distance from a medical device
including all connection cables**

Mobile phone	3.3 m
DECT telephone (cordless telephone)	1.2 m
Bluetooth devices (laptops, mobile phones)	0.7 m
Remote controls	1.2 m
WLAN devices (e.g.: laptops, repeater, access point, print server)	2.3 m

Warning

The use of accessories and cables other than those listed in the Instruction for Use can lead to increased interference emission, reduced interference immunity of the device or impairment of the application.

Warning

DO NOT use the device in an environment where there is contamination, smoke, gases, flammable anaesthetics, cleaning agents or chemical vapours.

Warning

Close supervision is required when this device is used in the vicinity of children or disabled persons.

Warning

A change in altitude above sea level can affect the performance of the device. Consult your physician or contact the manufacturer.

Warning

Stow and lay electrical cables properly to avoid potential tripping hazards.

Warning

Risk of infection if the device is used again!

If the device is used by several patients, infections can be transmitted to the next patient.

When the device is used again, it must be cleaned and disinfected, see chapter 5.

Warning

The ME equipment has no defibrillation protection.

Warning

The sensor is not protected against defibrillator discharge.

Warning

When using the device together with high-frequency surgical devices, the measuring sensor must be removed from the patient during a procedure with the HF surgical device in order to avoid burns to the patient and damage to the sensor.

Warning

Modification of the ME equipment is not permitted.

Warning

The device may not be modified without the manufacturer's permission.

Warning

If the ME equipment is modified, suitable inspections and tests must be carried out to ensure continued safe use, suitable inspections and tests must be carried out.

Warning

The device is not a blood gas analyser.

Warning

The sensor surface must not be covered with adhesive tape or similar.

Warning

Do not expose the sensors to direct sunlight or UV radiation.

Warning

Do not insert or remove the sensor modules except ONLY in the de-energized state.

Warning

Please avoid disconnecting the device from the computer while the live transmission is running. It is recommended to stop the transmission first and then remove the USB cable.

Warning

The minimum distance between the patient and the computer must be at least 1.5 m.

Warning

If the mains supply is interrupted, the ME equipment has no power supply.

Warning

The device must not be operated in a pressure chamber. If a measurement is carried out in a pressure chamber, only the sensor and the accessories, such as the interface box, may be in the pressure chamber.

2.6 Expected service life

If maintenance is carried out to the prescribed extent and at the specified intervals, the device will operate safely during the intervening period.

With each maintenance and Technical Safety Control, the "expected service life" according to IEC 60601-1 will therefore be prolonged until the next prescribed Technical Safety Check.

2.7 Disclaimer of liability

The device has been approved for use with certain accessories and spare parts (see Chapter 10).

The manufacturer accepts no responsibility or liability for personal injury or other damage and excludes any warranty for damage to the device resulting from the use of non-approved or unsuitable accessories or spare parts

2.8 Guarantee / Warranty

medicap homecare GmbH provides a two-year guarantee from the date of delivery for defects caused by material or manufacturing faults.

Defects covered by the warranty will be rectified in accordance with the warranty conditions.

In addition, medicap homecare GmbH does not provide a guarantee if the operator jeopardizes the functions of the device by failing to observe these Instructions for Use, by using the device for purposes other than those for which it is intended or by tampering with it.

medicap homecare GmbH does not automatically consider the owner of the device to be authorized to carry out service work.

Note

The guarantee can only be claimed in conjunction with the proof of purchase.

2.9 Adresses

Manufacturer:

medicap homecare GmbH
Hoherodskopfstr. 22
35327 Ulrichstein
Germany

Tel.: +496645/970-0
Fax: +496645/970-200

CE
0044

Service:

3 Operation

3.1 Preparation

3.1.1 Connecting the power supply unit

The power supply socket is located on the right side of the housing.

Plug the power supply into the contact socket (100 to 240 VAC).

Caution

Only the power supply unit supplied by medicap homecare GmbH may be used.

Caution

The built-in rechargeable battery enables a runtime of approx. 3 hours without any external power supply.

Caution

The battery must be protected against deep discharge. There is a specially designed deep discharge protection for this purpose. This automatically switches off the load as soon as the battery voltage falls below a set limit value. As the battery should not be more than 80% discharged, charging the battery in a timely manner prevents a deep discharge.

Caution

With the interface box, the device is also suitable for measurement in a pressure chamber. The device must NOT be in a pressure chamber, only the interface box.


Warning

The device must not be operated in a pressure chamber. If a measurement is carried out in a pressure chamber, only the sensor and the accessories, such as the interface box, may be in the pressure chamber.

3.1.2 Assignment of the sensors

Up to eight sensors can be connected to the Précise 8008 at the same time. Each sensor is identified by color and number. The respective color and number can be found in the display of the device.

The assignment of the sensors is shown visually in the following illustration:

Caution

If a sensor is activated but the sensor module is not connected, a measurement cannot be started!

3.1.3 Connecting sensors

Caution

The sensor modules may only be inserted or removed when they are de-energized.

Caution

The sensors must not be exposed to direct sunlight or UV radiation.

The required number of sensors can be connected laterally on the left and right side. A maximum of four sensors can be connected on each side. A maximum of 8 sensors can be connected to the Précise 8008.

If the sensor is activated in the device but not plugged in, this is indicated by ****. Sensors can also be deactivated manually. Deactivated sensors are indicated by "hyphens" in the sensor field.

Asterisk:
Sensor is activated,
but no sensor
module is
connected

Hyphens:
Sensor is
deactivated

Caution

If a sensor is activated but the sensor module is not connected, a measurement cannot be started!

3.1.4 Fixation & use of the supporting frame with a retaining arm


Placement of the 180° double-jointed arms into the Précise 8008 base unit:

The 180° double-jointed arms are fixed to the base unit of the device with a click system.

Carefully insert the optical fibre cable into the guide clips.

Movement options of the double-jointed 180° pivotable cable carrier:

3.2 Operation

3.2.1 General

Touching the screen: To select various menu items.

Touching the symbols: When switched on, the software is called up or confirmed by touching the symbols.

Back:  To exit the current menu without changes

Discard:  Directly back to the initial display, current data or selection is discarded

System settings:  System menu, for setting the language and system settings (service menu)


Note

If --- (hyphens) are displayed instead of the temperature, the sensor is deactivated.

3.2.2 Switching on the device

The device must be switched on with the I/O button (at least 3 sec.). This is located on the right-hand side of the device (acoustic signal sounds).

After switching on, the **initial screen** appears. The Précise 8008 is ready for operation when the tcpO₂ value of the optical sensors has risen to the air tcpO₂ value of the built-in control barometer.

3.2.3 Fixing the sensor and preparing the measurement site


After switching on the device, a short introduction on how to fix the sensor appears on the initial display.

This must then be confirmed with  in order to continue working.

The sensor can be fixed to various body parts, such as arms, legs, chest area, etc. The fixation site should be determined by a physician depending on the area of application.


The sensor surface should be cleaned before each measurement with a disinfectant¹ or alcohol swab². Do not exert too much pressure.


¹ Recommended disinfectants, see chapter 5.3 Cleaning and Disinfection, e.g. DESCOTON FORTE (conc. max. 4%), or SEKUSEPT AKTIV (Conc. max. 20g/l)

² Recommended disinfectants, see chapter 5.3 Cleaning and Disinfection, e.g. Alcohol swab from B.Braun

If necessary, clean the area of skin to be measured with an alcohol swab. Furthermore, the desired skin area must be depilated and dead epithelial cells removed by peeling off the skin several times with a new adhesive strip each time.

Carefully remove one adhesive ring from the backing material. The adhesive surface must not be soiled in the process.

Apply the adhesive ring to the skin area to be measured. Carefully drip a drop of contact liquid into the opening of the adhesive ring. Make sure that the adhesive ring is airtight and firmly attached to the skin.


The sensor head can be inserted into the adhesive ring using a one-touch click system. To prevent the sensor from shearing off, ensure that there is no tension on the cable when inserting it.

3.2.4 General operation - Patient menu and patient database

Touch the symbol for the Patient menu in the *initial display*:

The display **Patient menu** appears. Select the required menu item and confirm by touching the symbol.

Create patient data

Change patient data


Open patient database


Exit / Back


Create **patient data display:**


To create a patient, the mandatory fields first name , family name  and date of birth  must be filled in. The selection of gender  is optional (default preselection is male).

Entering the mandatory fields unlocks the „Create patient data“ button. The symbol  must be touched to create and save the patient in the database.


Change patient data display:

To change patient data, the mandatory fields first name , family name , date of birth  must be completed. All entered data can be changed.

Entering the mandatory fields unlocks the "Change patient data" button. The symbol  must be touched. The changes entered are now saved.

Display patient database:

Select the desired patient and confirm by touching the  symbol.

Note

By entering the first name or family name, you can quickly and easily narrow down the selection and save yourself the trouble of manually searching the database.

3.2.5 Measurement


The measurement graph, temperature and oxygen partial pressure can be read on the ***initial display***.

Record measurement display:

Oxygen partial pressure

Measurement duration

Current patient


The measurement is started by touching the symbol

Once the measurement has started, the symbol changes to

Display during the measurement:

Stability value detection

Stability value detection can be selected manually by touching the symbol.

As long as the bar lights up orange, the stability value detection is within the defined range until a stability value is reached.


Stop

The measurement can be stopped at any time by touching the symbol.

When the stability value is reached, an acoustic signal sounds and a line is displayed visually.

When the stability value is reached, the symbol on the display changes as follows:

The stability value detection process can be repeated a maximum of five times for provocations or position changes. An „event“ is always restarted by touching the symbol .

When an "event" is activated, the time of the stability value detection is displayed in color.

Scaling

By selecting the scaling between 3, 6 or 9 minutes, the last minutes of the measurement graph of a live measurement can be shown on the display.

If you touch the symbol  , the entire measurement is displayed again

End measurement

The measurement can be stopped at any time by touching the symbol  . Measurement data is automatically saved in the database.

All measurement data and measurement parameters are saved. The data can be called up as often as required.

Display Measurement has been completed:

Note

Once the measurement has been completed, scaling is no longer possible.

3.2.6 Sensor removing

After the measurement, the sensors must be removed from the skin. The **initial display** shows a brief introduction on how to remove and clean the sensor.

This must then be confirmed with  in order to continue working.

Remove the sensor from the adhesive ring and carefully pull the adhesive ring off the skin. The sensor must be cleaned with a disinfectant¹ or alcohol swab² after each measurement. Do not exert too much pressure.

¹ Recommended disinfectants, see chapter 5.3 Cleaning and Disinfection, e.g. DESCOTON FORTE (conc. max. 4%), or SEKUSEPT AKTIV (Conc. max. 20g/l)

² Recommended disinfectants, see chapter 5.3 Cleaning and Disinfection, e.g. Alcohol swab from B.Braun


3.2.7 Evaluate measurement data

If a measurement has just been completed, you can touch the symbol


to go directly to the database of the current patient.

To evaluate the measurement data of a patient from the database, you must first open the patient menu. Touch

the patient symbol  in the ***initial display*** to access the "Patient menu" selection menu. The patient must be selected.

Select the desired measurement and confirm by touching the symbol . The measurements are saved by date, latest first and time. The measurement graph and the associated data are called up and displayed graphically.

Note

If a USB stick is connected, the symbol appears on the display. If a printer is connected, the symbol appears on the display.

3.2.8 System menu

3.2.8.1 System settings

System settings display:

Active event list	
Software	v3.2 (11.05.2022)
DB-Text	v2.0.0 (28.03.2024)
DB-Measurement data:	v2.0 (07.02.2020)
Rechargeable battery status	97 %
Rechargeable battery voltage	15.6 V
<hr/>	
Controller	Feb 21 2020
Sensor 1	v0.13
Sensor 2	v0.13
Sensor 3	v-..
Sensor 4	v0.13
Sensor 5	v-..
Sensor 6	v-..
Sensor 7	v-..
Sensor 8	v-..
<hr/>	
Device-ID	808.010.24.0100
Device-IP	
<hr/>	
Default printer	DeskJet-3700-series
<hr/>	
Language	

System information

Delete data

Sound ON

Sound OFF

Back / Exit

Printer settings

Symbol printer functions

The printer function is described in more detail in section 3.3.

System settings display:

Active event list	
Software	v3.2 (11.05.2022)
DB-Text	v2.0.0 (28.03.2024)
DB-Measurement data:	v2.0 (07.02.2020)
Rechargeable battery status	97 %
Rechargeable battery voltage	15.6 V
<hr/>	
Controller	Feb 21 2020
Sensor 1	v0.13
Sensor 2	v0.13
Sensor 3	v-..
Sensor 4	v0.13
Sensor 5	v-..
Sensor 6	v-..
Sensor 7	v-..
Sensor 8	v-..
<hr/>	
Device-ID	808.010.24.0100
Device-IP	
<hr/>	
Default printer	DeskJet-3700-series
<hr/>	
Language	                    <img alt="Flag


The current software version, hardware parameters and serial number of the sensors can be read in this menu.

This menu is also used for the language setting, which can be selected using the drop-down menu. The required language can be selected in the form of the country flag.

After selecting the language, the device restarts automatically. The device does not need to be switched off

again. Once the language has been selected, exit the settings menu via the "Exit / Back" symbol .

Next, a new window will appear with the message that the software is being restarted (in the newly selected language). Shortly afterwards, the software appears in the new language.


Note

The user may only use the device in the language in which the Instruction for Use is provided.

3.2.8.2 Delete data

Display Delete data:

The interface consists of a vertical sidebar with icons and a main content area.

Top Screenshot: Shows a list of patients. A row for "John Doe" (12.09.1980) is selected, indicated by a blue border. To the right of the list are three icons: a trash can (Delete patient), a speaker (Cancel selection), and a gear (Select patient).

Bottom Screenshot: Shows a detailed view of measurements for the selected patient. The measurement table has columns for Date, Time, and a checkbox. The first measurement (02.04.2024, 14:38:59) has its checkbox checked. To the right of the table are two sets of icons: a trash can (Delete measurement) and a speaker (Cancel selection).

Annotations:

- Select patient:** Points to the gear icon in the top screenshot.
- Patient is selected:** Points to the checked checkbox in the bottom screenshot.
- Delete patient:** Points to the trash can icon in the top screenshot.
- Cancel selection:** Points to the speaker icon in the top screenshot.
- Select measurement:** Points to the gear icon in the bottom screenshot.
- Delete measurement:** Points to the trash can icon in the bottom screenshot.
- Cancel selection:** Points to the speaker icon in the bottom screenshot.

In the menu, individual **measurements of a patient or the patient itself** can be deleted.

Caution

The data cannot be restored after deletion.

To delete the data, click on the patient or the respective measurement in the settings menu and confirm with

. A window opens to confirm the deletion again:

Really delete the patient?

Really delete measurement data?

Touch the “Confirm“ button to permanently delete the selected data and return to the initial screen.

Touch the “Cancel“ button to return to the initial screen, without changing anything.

Data has been deleted!


Note

If a patient is deleted, all their measurements are also permanently deleted.

3.2.9 Date & time

The Select Menu for the date and time can be opened by touching the symbols on the *initial display*.

Touch the desired field to open a keyboard. The entry can be deleted and re-entered using the symbol .



Touch the symbol to confirm the entry. The symbol exits the entry and does not change it.

Touch the “Confirm“ button to save the changes you have made and return to the initial screen.

Touch the “Cancel“ button to return to the initial screen, without changing anything.


Note

The device restarts, so the process takes a little time. Please be patient.

3.2.10 Switching off the device

Switch off by touching the I/O button. This is located on the right-hand side of the device (acoustic signal sounds).


After confirming the query, the Précise 8008 is switched off.

3.3 Print functions

3.3.1 USB interface

The Précise 8008 has two USB interfaces. The USB-A and USB-B interface.

With the USB-B interface, the device offers the option of transferring the live measurement data to the associated visualization software, where it can be displayed and saved. In addition to the above-mentioned storage in the device, this variant offers the option of further processing the data.

The USB-A interface offers the option of exporting previously saved measurements to a USB stick for further external processing. The exported data can also be subsequently imported into the visualization software for further processing.

Note

The two usage variants are also presented in section 1.4 Brief description.

Note

A USB hub is required to connect a printer and a USB stick to the device at the same time. This is plugged into the USB-A port. The USB stick and the printer are plugged into the free slots of the USB hub.

3.3.2 Define settings


Touch the "Settings" button  in the main interface.

The settings menu opens. Click on the "Set printer settings" button: 

A submenu opens. The available printers are determined (including network printers). Select the appropriate printer and confirm it with the "Accept" symbol .

The currently selected printer is displayed above the list. You can exit the setting again using the "Exit / Back" symbol .

3.3.3 Print measurement

Once the above settings have been made, a new symbol appears on the main interface after a recording or a loaded measurement.

→ Touch Symbol "Print PDF"

If only one printer is connected, this printer appears on the main interface. If this is clicked, a PDF is generated in the background and the selected printer prints the document.

If only one USB stick is inserted, the "Save measurement as PDF on USB stick" symbol appears

After touching the symbol, a PDF is generated in the background and then saved to the connected USB stick.

A window opens after clicking on the respective symbol. Depending on the variant, the corresponding symbol lights up. If both a printer and a USB stick are inserted, both operations are carried out. This process may take some time. Once both processes have been completed, the window closes again.

In this case, only one printer is connected to the device.

3.3.4 Data export

The "Data export" symbol  is active when a USB stick is connected. If it is selected, the following display appears:

Symbol „Note“ not active

This can be used to save data to the USB stick. To do this, touch the "Measurement data" symbol  to open it.

To export measurements and measurement data from the device and import and display them in the appropriate visualization software, for example, select the patient from the first list.

Note

The export function is only optional. If the Précise 8008 device and the visualization software are used, the patient data and the measurement data are saved both in the device and in the visualization software. It is no longer necessary to export the data from the Précise 8008 device.


Note


The measurement can also be started in the visualization software, in which case the patient data and the measurement data are automatically saved in the visualization software.

After selection, all measurements available for the selected patient are displayed in the second list.

A measurement can be selected and the data exported by touching the  symbol.

As soon as the saving process is complete, the following display appears (display with a sample file):


This must be confirmed by touching the symbol .


Note

The exported file contains the patient data (first name, family name, date of birth; each encrypted; gender), the measurement parameter data (measurement date, measurement time) and the software version number of the device. The file also contains the times at which the stability values and events were set and the corresponding event names (default: E1, E2, E3, E4, E5). At the end, all measured values (max. 8 sensors) are saved in the file. With the help of this file, the measurement can be reconstructed in the visualization software / PC software and displayed again.

If this symbol  is active, there is an error report. This can be transferred to the connected USB stick. This is done by clicking on the symbol . The following display then appears:

To complete the process, the display must be confirmed by touching the symbol .

3.4 Sensor menu

Select the required sensor in the ***initial display*** and touch the "Sensor" symbol to open the "Sensor menu" (For example: sensor 2: ). Touch  to confirm the changes you have made and return to the initial screen.

To exit the menu without making any changes, use .

3.4.1 Set sensor temperature

In the "Parameters" menu, the Target temperature can be changed by touching-the  and  symbols.

Target temperature: required temperature of the sensor
 Actual temperature: current temperature of the sensor

Temperature range from 40°C - max. 44°C

Note

Reducing the temperature extends the measurement duration.

3.4.2 Switching the sensor on/off

The desired sensor can be selected in the ***initial display***. You can use the on/off symbol to activate or deactivate the sensor.

Activate (ON)

Deactivated (OFF)

The deactivated sensors are indicated by "hyphens" in the ***initial display***. This setting is not possible during an ongoing measurement.

Caution

If a sensor is activated but the sensor module is not connected, a measurement cannot be started!

If a sensor (e.g. with sensor 3) is set to active but not recognized (e.g. because it is not plugged in correctly), the following field message appears:

4 Malfunctions and their elimination

The Précise 8008 contains a microcontroller that ensures permanent monitoring of the most important parameters.

If the sensor temperature rises above 45°C, the heating is switched off and the display shows ::: (colons) instead of the temperature.

The device must be switched off. After the sensor has cooled down, the device must be switched on again. If the error occurs again, the technical service must be informed.

If the drift of +/- 6 mmHg is exceeded, recalibration is recommended.

Audible alarm signal (beep) indicates that a sensor is not connected.

Caution

If there is moisture, splash water or other contamination on the display, the display can no longer be operated.

Caution

If a sensor is activated but the sensor module is not connected, a measurement cannot be started!

5 Cleaning Process

Equipment cleaning is a process that a used product undergoes so that it can be safely reused. This procedure includes cleaning and disinfection.

In addition, a test and restoration of the technical and functional safety of the product used must be carried out.

Warning

Risk of infection if the device is used again! If the device is used by several patients, infections can be transmitted to the next patient. If the device is used again, it must be cleaned and disinfected.

5.1 General Notes

Note

Wear suitable protective equipment during cleaning / disinfection.

Caution

The cleaning and disinfecting agent must be used in accordance with the manufacturer's instructions, see the manufacturer's Instructions for Use.

Caution

Do not allow any liquid to enter the interior of the device.

Caution

DO NOT allow detergent to drip into the air inlet or outlet ports.

Caution

After cleaning, disinfection and a functional check, the device is suitable for reuse with further patients.

Warning

Risk of injury due to electric shock! Penetrating liquids from cleaning agents and disinfectants can cause a short circuit, injure the user and damage the device.

Caution

Ensure spatial separation of equipment that is still contaminated and equipment that has already been disinfected (such as work surfaces and storage rooms, etc.).

Caution

Switch off the device and disconnect it from the power supply before starting cleaning.

Caution

Do not open the device! Repairs may only be carried out by qualified technicians.

Caution

Do not pour liquid over the device and components or immerse them in liquids.

Caution

DO NOT clean the housing with concentrated chlorine-based products (ethylene chloride), oil-based products (Pine-Sol®, Lestoil®), or other aggressive chemical cleaning agents.

Use commercially available cleaning agents to clean the outer housing, but do not use abrasive, scratching or corrosive agents.

Caution

No lubricants other than those recommended by medicap homecare GmbH may be used.

5.2 Cleaning Intervals

Intervals	Action
Daily	- Clean / disinfect sensor
Monthly	- Clean / disinfect housing - Clean / disinfect display
Annual	- the device must be checked annually by personnel authorized by medicap homecare GmbH as part of the Technical Safety Control. - the sensors must be checked annually by personnel authorized by medicap homecare GmbH as part of the Technical Safety Control.
If required	- Clean / disinfect housing - Clean / disinfect display - Clean / disinfect sensor - Replace / Change approved disposables
For patient change	- Clean / disinfect housing - Clean / disinfect display - Clean / disinfect sensor - Replace / Change approved disposables

5.3 Cleaning and Disinfection

Component / part	Cleaning (manual)	Disinfection (manual)
Housing incl. handle	Wipe with a damp cloth; Use water or commercially available cleaning agents, but do not use abrasive, scratching or corrosive agents, allow to dry.	Wipe disinfection (see recommended disinfectants below), allow to dry.
Display (Touchscreen)	Spread the cleaning agent evenly on a clean, dry cloth. You can now use this cloth to clean the touchscreen surface. First clean the dirt from the edge of the surface towards the centre of the monitor and then remove the final dirt in the centre itself. Do not use abrasive, scratching or corrosive agents, allow to dry.	Spread the disinfectant evenly on a clean, dry cloth. You can now use this cloth to clean the touchscreen surface. First clean the dirt from the edge of the surface towards the centre of the monitor and then remove the final dirt in the centre itself. Do not use abrasive, scratching or corrosive agents (see recommended disinfectants below), allow to dry.
tcpO2 sensor	Do not use any pointed or sharp objects! The sensor surface must not be covered with adhesive tape or the like. The sensor must be cleaned with an alcohol swab after each measurement. Do not exert any great pressure on the white sensor surface! Allow the sensor to air dry.	Use alcohol swabs! Do not exert any significant pressure on the white sensor surface! Allow the sensor to air dry.
Power supply unit / power cable	Wipe with a damp cloth; Use water or commercially available cleaning agents, allow to dry.	Wipe disinfection (see recommended disinfectants below), allow to dry.
Adhesive rings	Single-use items are to be disposed of, no cleaning.	Single-use items are to be disposed of, no disinfection.

Cleaning of the devices is always done manually. The cleaning agent is to be applied on a cloth (cloths with low particle release are to be used).

Brushes, spray guns or similar are not recommended.

Caution

If there is moisture, splash water or other contamination on the display, the display can no longer be operated.

Recommendation for cleaning agents and disinfectants

The following disinfectants have been validated by medicap homecare GmbH and can be used:

- Biguanid Fläche N (Dr. Schuhmacher)
- Descogen Liquid
- microzid AF wipes (Schülke)
- KORSOLEX FF Konzentrat

Recommendation for cleaning agents and disinfectants for the tcpO2 sensor:

- DESCOTON FORTE (Conz. max. 4%)
- SEKUSEPT AKTIV (Conz. max. 20g/1Liter)
- Alcohol swab B.Braun REF 9160612 / PZN 00629703
- Mikrozid ® sensitive wipes (Schülke)

Any other disinfectant used must be listed in the VAH list and be compatible with the materials to be treated.

Drying is carried out in fresh air without certain drying media (such as hair dryers). Wiping with cloths (for drying) is not recommended.

6 Maintenance

Caution

All maintenance intervals are guide values. In individual cases, these may be shorter due to changes in environmental conditions with regard to temperature, humidity, dust, etc.

Caution

Maintenance may only be carried out by an authorized service technician.

Caution

Circuit diagrams, component lists, descriptions or other information can be provided on request.

Caution

During maintenance, only monitored test equipment and measuring devices suitable for the intended purpose may be used.

The device and the sensors must be checked and serviced annually by the supplier or manufacturer or by authorized trained personnel.

For commercially used devices, a safety check must be carried out.

The device must not be maintained or serviced during use on the patient!

Caution

Works on the opened device may only be carried out by an authorized service technician.

7 Pictorial Symbols

7.1 Definition of the symbols on the device or accessories

Attention! Caution!

Follow the Instructions for Use

Fragile

Keep away from sunlight

Store in a dry place

Follow the Instructions for Use

protection class II

Single-use; Do not reuse

Article number

Batch number

Serial number

The device must not be disposed of with household waste

CE- Mark 0044; Notified body: TÜV NORD CERT GmbH 0044

Degree of protection against electric shock: Type BF

Alternating current

Dangerous electrical voltage

Manufacturer

Date of manufacture

Medical device



UDI number - Unique device identifier

European Authorized Representative

7.2 Software symbols

	Patient menu		Wiping away the data / cancel selection
	Create patient data		Measurement data
	Select patient		Date
	Patient database		Database
	Settings		Event - Set measurement point
	Date of birth		Gender
	Note / Caution		Info
	Accept		Cancel
	First name		Family name
	Record measurement		Clean sensor

	Stability value detection		Temperature
	Start measurement		Stop measurement
	Data export - save data to USB stick		Eject USB port
	Exit / back		Barometer
	Delete		Export data
	Sound ON		Sound OFF
	Measurement duration		Setting the scale
	Rechargeable battery empty (under 5%)		Rechargeable battery very low (under 15%)
	Rechargeable battery low (under 25%)		Rechargeable battery low (under 25%) in charge level
	Rechargeable battery defective		Rechargeable battery fully charged
	Rechargeable battery empty in charge level		Battery 50% charge level
	Save measurement as PDF on USB stick		Define printer settings
	Waiting for process conclusion		Print measurement as PDF file
	Signals a connection to the visualization software		Sensor number

Sensor 3 set to active,
temperature sensor
defective

Sensor 3 set to active but
not detected

8 Technical Data and Specifications

Functionality is guaranteed at room temperature, normal humidity (50 - 55 %), good lighting conditions (at least 750 lux) and no background noise.

The ambient volume should not exceed room volume, up to a maximum of 55 dB.

Dimensions:	approx. 400 mm x 250 mm x 170 mm
Weight (inkl. 8 Sensors):	approx. 3920 gr.
Display:	15,1"
Operating voltage:	100 to 230 VAC +/- 10%,
max. Power consumption:	35 VA
Protection class:	II
Degree of protection against electric shock:	Type BF
MDD 93/42 EEC classification:	IIa
Device corresponds:	MDD 93/42/EWG
Ambient temperature Operating:	15 bis 35°C
Relative humidity:	non-condensing 10 to 95 %
Storage temperature:	-10 to 50°C
Oxygen partial pressure:	0 to 2000 mmHg +/- 10% mmHg
Adjustable sensor temperature:	40 to 44°C
No cross-sensitivity:	pH 1 – 14, CO2, H2S, SO2, Ionic species
Cross-sensitivity:	Organic solvents such as acetone, toluene, chloroform or methylene chloride, chlorine gas
Degree of protection against penetration of liquids:	IP 22
Notified body:	 CE 0044, TÜV Nord Cert

9 Storage and Scope of delivery

The device is delivered safely packed in the box. Transport and storage take place in the original box. Store the device, before and after use, under the prescribed environmental conditions. If the device has been used, the device should be cleaned and disinfected before storage.

10 Recommended Accessories, Disposables and Spare Parts

The accessories, disposables or spare parts listed below have been validated by medicap homecare GmbH for use with the Précise 8008 medical device and are recommended.

Caution

The accessories such as the adhesive rings are "single use" items. They are NOT suitable for repeated use by a single patient as well as by multiple patients in the physician's office/clinic.

Caution

The use of accessories, spare parts or disposables that have NOT been approved may result in impaired performance.

Caution

The specified disposables have been validated by medicap homecare GmbH for use with the devices. If disposables from other manufacturers are used, the measurement accuracy can no longer be guaranteed.

Caution

The specified printer has been validated by medicap homecare GmbH for use with the device. If printers from other manufacturers are used, the printing performance or printer accuracy can no longer be guaranteed.

Article or spare part number	Description
802.204	Adhesive rings for sensor (48 pcs.)
802.280	Interface-Box
802.218	Printer for P8008 HP DeskJet 3762
802.203	Contact liquid 15ml bottle
808.060	tcpO2-Sensor P800x Nr.1 RED
808.061	tcpO2-Sensor P800x Nr.2 YELLOW
808.062	tcpO2-Sensor P800x Nr.3 ORANGE
808.063	tcpO2-Sensor P800x Nr.4 GREEN
808.064	tcpO2-Sensor P800x Nr.5 BROWN
808.065	tcpO2-Sensor P800x Nr.6 GRAY
808.066	tcpO2-Sensor P800x Nr.7 BLACK
808.067	tcpO2-Sensor P800x Nr.8 LILA/ PURPLE
801.312	Power supply Précise 8001 / Précise 8008 (Gen.2) 15V + Länderstecker
802.040	Précise 800X Visualization Software Pack 2.x

11 Disposal

medicap homecare GmbH complies with the criteria of the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU.

The product contains substances that may harm the environment if not disposed of in accordance with national legislation.

Therefore, the device or packaging can be returned to medicap homecare GmbH for disposal free of charge. We ensure environmentally friendly disposal.

Do not dispose of used batteries with household waste!

Accessories that are not part of the product must be disposed of in accordance with the relevant disposal labelling.

medicap homecare GmbH
Hoherodskopfstr. 22
35327 Ulrichstein
Germany
Tel.: +49 6645/970-0
Fax: +49 6645/970-200

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
The screenshot shows a software interface for managing patient data. At the top, there are fields for 'nombre' (name) and 'apellido' (surname), each with a small blue circular icon. Below these are two tables:

nombre	apellido
Fulano	de Tal

hora	fecha
08:33:20	18.09.2024
14:59:07	17.09.2024
11:00:58	17.09.2024

At the top right, there are three icons: a blue circle with a white document and a gear, a blue circle with a white floppy disk, and a blue circle with a white document and a downward arrow. The entire interface has a light gray background.

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
The screenshot shows a software interface for managing patient data and measurements. At the top, there are two input fields labeled 'Nome' (Name) and 'Cognome' (Surname), both containing the placeholder text 'Nome' and 'Cognome'. Below these fields is a table with two columns: 'Nome' (Name) and 'Cognome' (Surname). The first row contains the values 'Mario' and 'Rossi'. To the right of this table is another table with two columns: 'Ora' (Time) and 'Data' (Date). It contains four rows of data:

Ora	Data
13:06:25	25.09.2024
13:04:25	25.09.2024
12:53:18	25.09.2024

At the top right of the interface, there are three circular icons representing different functions: a blue icon with a white downward arrow, a blue icon with a white floppy disk, and a blue icon with a white exit symbol.

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